

Dionex_COMPASS_TEMPEST

Stephanie J. Wilson

2023-06-09

Daily Set up

Read in and Format the raw data - change wd & file names (Make sure files are UTF-8 encoded)

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_TEMPEST_2023_S04a.txt", sep='\t' , header=T, skip=3)
head(Sdat)
```

```
##      X          X.1                X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank          Unknown    n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 5.733  0.4858  6.82  0.6526  3.36
## 6 6 Standard 2 Calibration Standard 5.767  0.9999  6.94  1.3433  6.78
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[,c(2,5,7)] # dont need this here
head(Sdat)
```

```
##          X.1 IC.S04.1 IC.S04.3
## 1 Lab Blank    n.a.    n.a.
## 2 Lab Blank    n.a.    n.a.
## 3 Lab Blank    n.a.    n.a.
## 4 Lab Blank    n.a.    n.a.
## 5 Standard 1  0.4858  0.6526
## 6 Standard 2  0.9999  1.3433
```

```
## Name the columns correctly
```

```
colnames(Sdat) <- c("Sample_ID", "S04_ppm", "S04_Area")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat$S04_Area <- as.numeric(Sdat$S04_Area)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm S04_Area
## 1 Lab Blank      NA      NA
## 2 Lab Blank      NA      NA
## 3 Lab Blank      NA      NA
## 4 Lab Blank      NA      NA
## 5 Standard 1  0.4858  0.6526
## 6 Standard 2  0.9999  1.3433
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("Raw Data/COMPASS_TEMPEST_2023_Cla.txt", sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 2 2 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 5 5 Standard 1 Calibration Standard 4.523  4.9712  93.18  8.9150  54.84
## 6 6 Standard 2 Calibration Standard 4.497  10.0411  93.01  18.0069  109.84
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[,c(2,5,7)]
head(Cldat)
```

```
##      X.1 IC.Cl.1 IC.Cl.3
## 1 Lab Blank  n.a.    n.a.
## 2 Lab Blank  n.a.    n.a.
## 3 Lab Blank  n.a.    n.a.
## 4 Lab Blank  n.a.    n.a.
## 5 Standard 1 4.9712  8.9150
## 6 Standard 2 10.0411 18.0069
```

```
## Name the columns correctly
```

```
colnames(Cldat) <- c("Sample_ID", "Cl_ppm", "Cl_Area")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat$Cl_Area <- as.numeric(Cldat$Cl_Area)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##      Sample_ID Cl_ppm Cl_Area
## 1 Lab Blank      NA      NA
## 2 Lab Blank      NA      NA
## 3 Lab Blank      NA      NA
## 4 Lab Blank      NA      NA
## 5 Standard 1  4.9712  8.9150
## 6 Standard 2 10.0411 18.0069
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1      154.3618  26.9769 1207.9926 359.0921
## 2      154.3618  26.9769  76.2651  22.4134
## 3      154.3618  26.9769   0.0009   0.0000
## 4      154.3618  26.9769 193.2321  69.7428
## 5      154.3618  26.9769      NA      NA
## 6      21.9166   2.1247 1207.9926 359.0921
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##              Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 26      100_TEMPEST_C_B4_20230810  3.8752  0.1041 0.0705  0.0025
## 27 100_TEMPEST_C_B4_20230810_dup  4.2038  0.1129 0.0396  0.0014
## 28      101_TEMPEST_C_H6_20230807  4.6586  0.1252 0.9717  0.0349
## 29 101_TEMPEST_C_H6_20230807_spk 91.2445  2.4515 1.0566  0.0379
## 30      102_TEMPEST_C_H3_20230807  5.4909  0.1475 0.0087  0.0003
## 31      103_TEMPEST_C_F4_20230810  4.8364  0.1299 0.0371  0.0013
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)
```

```
##              Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 26      100_TEMPEST_C_B4_20230810  3.8752  0.1041 0.0705  0.0025
## 27 100_TEMPEST_C_B4_20230810_dup  4.2038  0.1129 0.0396  0.0014
## 28      101_TEMPEST_C_H6_20230807  4.6586  0.1252 0.9717  0.0349
## 29 101_TEMPEST_C_H6_20230807_spk 91.2445  2.4515 1.0566  0.0379
## 30      102_TEMPEST_C_H3_20230807  5.4909  0.1475 0.0087  0.0003
## 31      103_TEMPEST_C_F4_20230810  4.8364  0.1299 0.0371  0.0013
```

```
all_dat[is.na(all_dat)] <- 0
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grep1("Standard", all_dat$Sample_ID),]
head(stds)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 297 Standard 1  0.4858  0.6526 4.9712  8.9150
## 298 Standard 1  0.4858  0.6526 5.0469  9.0507
## 299 Standard 1  0.4858  0.6526 5.0406  9.0394
## 300 Standard 1  0.4924  0.6615 4.9712  8.9150
## 301 Standard 1  0.4924  0.6615 5.0469  9.0507
## 302 Standard 1  0.4924  0.6615 5.0406  9.0394
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(S04_ppm), sd = sd(S04_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID   mean      sd    cv flag
##   <fct>      <dbl>   <dbl> <dbl> <chr>
## 1 Standard 1  0.490 0.00293 0.599 YES
## 2 Standard 2  1.00 0.00216 0.216 YES
## 3 Standard 3  2.01 0.00351 0.175 YES
## 4 Standard 4 10.1  0.0244  0.243 YES
## 5 Standard 5 20.0  0.0554  0.277 YES
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID   mean      sd    cv flag
##   <fct>      <dbl>   <dbl> <dbl> <chr>
## 1 Standard 1   5.02 0.0364 0.725 YES
## 2 Standard 2  10.1 0.0540 0.534 YES
## 3 Standard 3  20.1 0.0868 0.431 YES
## 4 Standard 4 101.  0.226 0.225 YES
## 5 Standard 5 200.  0.220 0.110 YES
```

```
lmS <- lm(stds$S04_Area ~ stds$S04_ppm)
S04_sum <- summary(lmS)
S04_Slope <- S04_sum$coefficients[2, 1]
S04_Int <- S04_sum$coefficients[1, 1]
```

```
lmCl <- lm(stds$Cl_Area ~ stds$Cl_ppm)
Cl_sum <- summary(lmCl)
```

```
Cl_Slope <- Cl_sum$coefficients[2, 1]
Cl_Int <- Cl_sum$coefficients[1, 1]
```

Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##              Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 26      100_TEMPEST_C_B4_20230810  3.8752  0.1041 0.0705  0.0025
## 27 100_TEMPEST_C_B4_20230810_dup  4.2038  0.1129 0.0396  0.0014
## 28      101_TEMPEST_C_H6_20230807  4.6586  0.1252 0.9717  0.0349
## 29 101_TEMPEST_C_H6_20230807_spk 91.2445  2.4515 1.0566  0.0379
## 30      102_TEMPEST_C_H3_20230807  5.4909  0.1475 0.0087  0.0003
## 31      103_TEMPEST_C_F4_20230810  4.8364  0.1299 0.0371  0.0013
```

```
# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06      #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)
```

```
##              Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 26      100_TEMPEST_C_B4_20230810  3.8752  0.1041 0.0705  0.0025 0.1208734
## 27 100_TEMPEST_C_B4_20230810_dup  4.2038  0.1129 0.0396  0.0014 0.1311229
## 28      101_TEMPEST_C_H6_20230807  4.6586  0.1252 0.9717  0.0349 0.1453088
## 29 101_TEMPEST_C_H6_20230807_spk 91.2445  2.4515 1.0566  0.0379 2.8460543
## 30      102_TEMPEST_C_H3_20230807  5.4909  0.1475 0.0087  0.0003 0.1712695
## 31      103_TEMPEST_C_F4_20230810  4.8364  0.1299 0.0371  0.0013 0.1508546
##              Cl_mM      salinity
## 26 0.0019887165 0.0001533935
## 27 0.0011170663 0.0000975572
## 28 0.0274104372 0.0017818619
## 29 0.0298053597 0.0019352762
## 30 0.0002454161 0.0000417209
## 31 0.0010465444 0.0000930397
```

Pull out dups and check with percent difference

```
#Show me the data that we have from the calculations
head(sampledat)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 26      100_TEMPEST_C_B4_20230810  3.8752  0.1041 0.0705  0.0025 0.1208734
## 27 100_TEMPEST_C_B4_20230810_dup  4.2038  0.1129 0.0396  0.0014 0.1311229
## 28      101_TEMPEST_C_H6_20230807  4.6586  0.1252 0.9717  0.0349 0.1453088
## 29 101_TEMPEST_C_H6_20230807_spk 91.2445  2.4515 1.0566  0.0379 2.8460543
## 30      102_TEMPEST_C_H3_20230807  5.4909  0.1475 0.0087  0.0003 0.1712695
## 31      103_TEMPEST_C_F4_20230810  4.8364  0.1299 0.0371  0.0013 0.1508546
##           Cl_mM      salinity
## 26 0.0019887165 0.0001533935
## 27 0.0011170663 0.0000975572
## 28 0.0274104372 0.0017818619
## 29 0.0298053597 0.0019352762
## 30 0.0002454161 0.0000417209
## 31 0.0010465444 0.0000930397
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "_dup")) #have to change this to match data
head(dups)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1 100_TEMPEST_C_B4_20230810_dup  4.2038  0.1129  0.0396  0.0014 0.1311229
## 2 110_TEMPEST_SW_C3_20231002_dup 125.8279  3.3807 983.5062 35.2749 3.9247629
## 3 120_TEMPEST_SW_B4_202311_dup  57.1233  1.5348 351.4739 12.6061 1.7817623
## 4 130_TEMPEST_FW_F6_20231208_dup  5.1334  0.1379  0.5484  0.0197 0.1601185
## 5 140_TEMPEST_C_H3_20231211_dup  5.4889  0.1475  0.5758  0.0207 0.1712071
## 6 90_TEMPEST_FW_D5_20230807_dup  5.1056  0.1372  0.0475  0.0017 0.1592514
##           Cl_mM      salinity
## 1 0.001117066 0.0000975572
## 2 27.743475317 1.7772217034
## 3 9.914637518 0.6351393373
## 4 0.015469676 0.0010169588
## 5 0.016242595 0.0010664706
## 6 0.001339915 0.0001118325
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "_dup")) %>%
  filter(!str_detect(Sample_ID, "_spk"))
head(sampledat2)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1 100_TEMPEST_C_B4_20230810  3.8752  0.1041 0.0705  0.0025 0.1208734
## 2 101_TEMPEST_C_H6_20230807  4.6586  0.1252 0.9717  0.0349 0.1453088
## 3 102_TEMPEST_C_H3_20230807  5.4909  0.1475 0.0087  0.0003 0.1712695
## 4 103_TEMPEST_C_F4_20230810  4.8364  0.1299 0.0371  0.0013 0.1508546
```

```
## 5 104_TEMPEST_C_C6_20230807 4.8123 0.1293 3.7102 0.1331 0.1501029
## 6 105_TEMPEST_SW_I5_20231002 114.6553 3.0805 281.4327 10.0940 3.5762726
##      Cl_mM      salinity
## 1 0.0019887165 0.0001533935
## 2 0.0274104372 0.0017818619
## 3 0.0002454161 0.0000417209
## 4 0.0010465444 0.0000930397
## 5 0.1046600846 0.0067303314
## 6 7.9388631876 0.5085748889
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[ , -c(2:5)]
colnames(dups) <- c('Sample_ID', 'SO4_mM_dup', 'Cl_mM_dup', 'salinity_dup')
head(dups)
```

```
##      Sample_ID SO4_mM_dup Cl_mM_dup salinity_dup
## 1 100_TEMPEST_C_B4_20230810 0.1311229 0.001117066 0.0000975572
## 2 110_TEMPEST_SW_C3_20231002 3.9247629 27.743475317 1.7772217034
## 3 120_TEMPEST_SW_B4_202311 1.7817623 9.914637518 0.6351393373
## 4 130_TEMPEST_FW_F6_20231208 0.1601185 0.015469676 0.0010169588
## 5 140_TEMPEST_C_H3_20231211 0.1712071 0.016242595 0.0010664706
## 6 90_TEMPEST_FW_D5_20230807 0.1592514 0.001339915 0.0001118325
```

```
#put it back together with the old data set and look for duplicates
QAdups <- merge(sampledat2, dups)
head(QAdups)
```

```
##      Sample_ID SO4_ppm SO4_Area Cl_ppm Cl_Area SO4_mM
## 1 100_TEMPEST_C_B4_20230810 3.8752 0.1041 0.0705 0.0025 0.1208734
## 2 110_TEMPEST_SW_C3_20231002 126.0340 3.3862 988.1749 35.4423 3.9311915
## 3 120_TEMPEST_SW_B4_202311 57.4551 1.5437 353.8215 12.6903 1.7921117
## 4 130_TEMPEST_FW_F6_20231208 5.2059 0.1399 0.5455 0.0196 0.1623799
## 5 140_TEMPEST_C_H3_20231211 5.5834 0.1500 0.3123 0.0112 0.1741547
## 6 90_TEMPEST_FW_D5_20230807 5.0450 0.1355 0.0494 0.0018 0.1573612
##      Cl_mM      salinity SO4_mM_dup Cl_mM_dup salinity_dup
## 1 0.001988717 0.0001533935 0.1311229 0.001117066 0.0000975572
## 2 27.875173484 1.7856580443 3.9247629 27.743475317 1.7772217034
## 3 9.980860367 0.6393814505 1.7817623 9.914637518 0.6351393373
## 4 0.015387870 0.0010117185 0.1601185 0.015469676 0.0010169588
## 5 0.008809591 0.0005903261 0.1712071 0.016242595 0.0010664706
## 6 0.001393512 0.0001152658 0.1592514 0.001339915 0.0001118325
```

```
QAdups$SO4_dups_chk <- ((abs(QAdups$SO4_mM-QAdups$SO4_mM_dup))/((QAdups$SO4_mM+QAdups$SO4_mM_dup)/2))*100
QAdups$SO4_dups_flag <- ifelse(QAdups$SO4_dups_chk <10, 'YES', 'NO, rerun')

QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')

head(QAdups)
```

```
##      Sample_ID SO4_ppm SO4_Area Cl_ppm Cl_Area SO4_mM
```


```
## 1 100_TEMPEST_C_B4_20230810 3.8752 0.1041 0.0705 0.0025 0.1208734
## 2 110_TEMPEST_SW_C3_20231002 126.0340 3.3862 988.1749 35.4423 3.9311915
## 3 120_TEMPEST_SW_B4_202311 57.4551 1.5437 353.8215 12.6903 1.7921117
## 4 130_TEMPEST_FW_F6_20231208 5.2059 0.1399 0.5455 0.0196 0.1623799
## 5 140_TEMPEST_C_H3_20231211 5.5834 0.1500 0.3123 0.0112 0.1741547
## 6 90_TEMPEST_FW_D5_20230807 5.0450 0.1355 0.0494 0.0018 0.1573612
##      Cl_mM      salinity S04_mM_dup      Cl_mM_dup salinity_dup S04_dups_chk
## 1 0.001988717 0.0001533935 0.1311229 0.001117066 0.0000975572 8.1346701
## 2 27.875173484 1.7856580443 3.9247629 27.743475317 1.7772217034 0.1636611
## 3 9.980860367 0.6393814505 1.7817623 9.914637518 0.6351393373 0.5791668
## 4 0.015387870 0.0010117185 0.1601185 0.015469676 0.0010169588 1.4024160
## 5 0.008809591 0.0005903261 0.1712071 0.016242595 0.0010664706 1.7069624
## 6 0.001393512 0.0001152658 0.1592514 0.001339915 0.0001118325 1.1940181
## S04_dups_flag Cl_dups_chk Cl_dups_flag
## 1 YES 56.1307902 NO, rerun
## 2 YES 0.4735756 YES
## 3 YES 0.6657069 YES
## 4 YES 0.5302130 YES
## 5 YES 59.3401644 NO, rerun
## 6 YES 3.9215686 YES
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
    color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```

```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows   Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         6 NO, rerun         2     6      100   33.3333
## 2  YES         6   YES           4     6      100   66.6667
```

Pull out dups and check with cv

```

#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)

```

```

##           Sample_ID  S04_ppm S04_Area  Cl_ppm Cl_Area  S04_mM
## 1 100_TEMPEST_C_B4_20230810   3.8752   0.1041   0.0705   0.0025 0.1208734
## 2 110_TEMPEST_SW_C3_20231002  126.0340   3.3862  988.1749  35.4423 3.9311915
## 3 120_TEMPEST_SW_B4_202311   57.4551   1.5437  353.8215  12.6903 1.7921117
## 4 130_TEMPEST_FW_F6_20231208   5.2059   0.1399   0.5455   0.0196 0.1623799
## 5 140_TEMPEST_C_H3_20231211   5.5834   0.1500   0.3123   0.0112 0.1741547
## 6 90_TEMPEST_FW_D5_20230807   5.0450   0.1355   0.0494   0.0018 0.1573612
##           Cl_mM  salinity S04_mM_dup  Cl_mM_dup salinity_dup S04_dups_chk
## 1 0.001988717 0.0001533935 0.1311229 0.001117066 0.0000975572 8.1346701
## 2 27.875173484 1.7856580443 3.9247629 27.743475317 1.7772217034 0.1636611
## 3 9.980860367 0.6393814505 1.7817623 9.914637518 0.6351393373 0.5791668
## 4 0.015387870 0.0010117185 0.1601185 0.015469676 0.0010169588 1.4024160
## 5 0.008809591 0.0005903261 0.1712071 0.016242595 0.0010664706 1.7069624
## 6 0.001393512 0.0001152658 0.1592514 0.001339915 0.0001118325 1.1940181
## S04_dups_flag Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1 YES 56.1307902 NO, rerun 5.7520804 YES
## 2 YES 0.4735756 YES 0.1157259 YES
## 3 YES 0.6657069 YES 0.4095327 YES
## 4 YES 0.5302130 YES 0.9916579 YES
## 5 YES 59.3401644 NO, rerun 1.2070047 YES
## 6 YES 3.9215686 YES 0.8442983 YES

```


```

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar,ncol=2, nrow=1)

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows   Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         6 NO, rerun         2      6      100    33.3333
## 2  YES         6   YES           4      6      100    66.6667
```

Pull out spikes and check with dionex output conc.

#Show me the data that we have from the calculations

```
head(sampledat)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area  S04_mM
## 26 100_TEMPEST_C_B4_20230810 3.8752 0.1041 0.0705 0.0025 0.1208734
## 27 100_TEMPEST_C_B4_20230810_dup 4.2038 0.1129 0.0396 0.0014 0.1311229
## 28 101_TEMPEST_C_H6_20230807 4.6586 0.1252 0.9717 0.0349 0.1453088
## 29 101_TEMPEST_C_H6_20230807_spk 91.2445 2.4515 1.0566 0.0379 2.8460543
## 30 102_TEMPEST_C_H3_20230807 5.4909 0.1475 0.0087 0.0003 0.1712695
## 31 103_TEMPEST_C_F4_20230810 4.8364 0.1299 0.0371 0.0013 0.1508546
##           Cl_mM      salinity
## 26 0.0019887165 0.0001533935
## 27 0.0011170663 0.0000975572
## 28 0.0274104372 0.0017818619
## 29 0.0298053597 0.0019352762
## 30 0.0002454161 0.0000417209
## 31 0.0010465444 0.0000930397
```

#pull out any rows that have "spk" in the SampleID column

```
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk"))      #have to change this to match data
head(spks)
```

```
##           Sample_ID S04_ppm S04_Area  Cl_ppm Cl_Area  S04_mM
## 1 101_TEMPEST_C_H6_20230807_spk 91.2445 2.4515 1.0566 0.0379 2.846054
## 2 111_TEMPEST_FW_I5_20231002_spk 89.9420 2.4165 0.3888 0.0139 2.805427
## 3 121_TEMPEST_SW_H6_202311_spk 106.4642 2.8604 141.0457 5.0588 3.320780
## 4 131_TEMPEST_FW_E3_20231208_spk 90.8868 2.4419 0.8722 0.0313 2.834897
## 5 141_TEMPEST_C_C6_20231208_spk 88.8596 2.3874 1.1352 0.0407 2.771666
## 6 91_TEMPEST_FW_B4_20230807_spk 93.5863 2.5144 0.0059 0.0002 2.919099
##           Cl_mM      salinity
## 1 0.0298053597 0.0019352762
## 2 0.0109675599 0.0007285616
## 3 3.9787221439 0.2548955799
## 4 0.0246036671 0.0016020654
## 5 0.0320225670 0.0020773064
## 6 0.0001664316 0.0000366613
```

#remove the dup from these IDs so we will have duplicate sample names

```
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[ ,-c(2,3,4 ,5, 7,8)]
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
head(spks)
```

```
##           Sample_ID S04_mM_spk
## 1 101_TEMPEST_C_H6_20230807 2.846054
## 2 111_TEMPEST_FW_I5_20231002 2.805427
## 3 121_TEMPEST_SW_H6_202311 3.320780
## 4 131_TEMPEST_FW_E3_20231208 2.834897
## 5 141_TEMPEST_C_C6_20231208 2.771666
## 6 91_TEMPEST_FW_B4_20230807 2.919099
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledats, spks)
head(QAspks)
```

```
##           Sample_ID S04_ppm S04_Area   Cl_ppm Cl_Area   S04_mM
## 1  101_TEMPEST_C_H6_20230807  4.6586   0.1252   0.9717  0.0349 0.1453088
## 2  111_TEMPEST_FW_I5_20231002  4.4138   0.1186   0.3448  0.0124 0.1376731
## 3   121_TEMPEST_SW_H6_202311 22.0449   0.5923 142.0156  5.0936 0.6876138
## 4  131_TEMPEST_FW_E3_20231208  7.6841   0.2065   0.8009  0.0287 0.2396787
## 5  141_TEMPEST_C_C6_20231208  6.1506   0.1653   0.8842  0.0317 0.1918465
## 6   91_TEMPEST_FW_B4_20230807  4.9651   0.1334  32.8252  1.1773 0.1548690
##           Cl_mM      salinity S04_mM_spk
## 1 0.027410437 0.0017818619   2.846054
## 2 0.009726375 0.0006490536   2.805427
## 3 4.006081805 0.2566481892   3.320780
## 4 0.022592384 0.0014732263   2.834897
## 5 0.024942172 0.0016237494   2.771666
## 6 0.925957687 0.0593411364   2.919099
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10              # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol      # mmoles of S04
head(QAspks)
```

```
##           Sample_ID S04_ppm S04_Area   Cl_ppm Cl_Area   S04_mM
## 1  101_TEMPEST_C_H6_20230807  4.6586   0.1252   0.9717  0.0349 0.1453088
## 2  111_TEMPEST_FW_I5_20231002  4.4138   0.1186   0.3448  0.0124 0.1376731
## 3   121_TEMPEST_SW_H6_202311 22.0449   0.5923 142.0156  5.0936 0.6876138
## 4  131_TEMPEST_FW_E3_20231208  7.6841   0.2065   0.8009  0.0287 0.2396787
## 5  141_TEMPEST_C_C6_20231208  6.1506   0.1653   0.8842  0.0317 0.1918465
## 6   91_TEMPEST_FW_B4_20230807  4.9651   0.1334  32.8252  1.1773 0.1548690
##           Cl_mM      salinity S04_mM_spk S04_spk_Conc
## 1 0.027410437 0.0017818619   2.846054 7.797879e-05
## 2 0.009726375 0.0006490536   2.805427 7.797879e-05
## 3 4.006081805 0.2566481892   3.320780 7.797879e-05
## 4 0.022592384 0.0014732263   2.834897 7.797879e-05
## 5 0.024942172 0.0016237494   2.771666 7.797879e-05
## 6 0.925957687 0.0593411364   2.919099 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
```

```

QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

```

```

##           Sample_ID S04_ppm S04_Area   Cl_ppm Cl_Area   S04_mM
## 1 101_TEMPEST_C_H6_20230807 4.6586 0.1252 0.9717 0.0349 0.1453088
## 2 111_TEMPEST_FW_I5_20231002 4.4138 0.1186 0.3448 0.0124 0.1376731
## 3 121_TEMPEST_SW_H6_202311 22.0449 0.5923 142.0156 5.0936 0.6876138
## 4 131_TEMPEST_FW_E3_20231208 7.6841 0.2065 0.8009 0.0287 0.2396787
## 5 141_TEMPEST_C_C6_20231208 6.1506 0.1653 0.8842 0.0317 0.1918465
## 6 91_TEMPEST_FW_B4_20230807 4.9651 0.1334 32.8252 1.1773 0.1548690
##           Cl_mM      salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 0.027410437 0.0017818619 2.846054 7.797879e-05      50 0.001483
## 2 0.009726375 0.0006490536 2.805427 7.797879e-05      50 0.001483
## 3 4.006081805 0.2566481892 3.320780 7.797879e-05      50 0.001483
## 4 0.022592384 0.0014732263 2.834897 7.797879e-05      50 0.001483
## 5 0.024942172 0.0016237494 2.771666 7.797879e-05      50 0.001483
## 6 0.925957687 0.0593411364 2.919099 7.797879e-05      50 0.001483

```

#gives us the total SO4 in the sample in mmoles

```
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)
```

##total SO4 in spiked sample in mmoles

```
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)
```

```
QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
```

```
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')
```


```
head(QAspks)
```

```
##           Sample_ID S04_ppm S04_Area   Cl_ppm Cl_Area   S04_mM
## 1 101_TEMPEST_C_H6_20230807 4.6586 0.1252 0.9717 0.0349 0.1453088
## 2 111_TEMPEST_FW_I5_20231002 4.4138 0.1186 0.3448 0.0124 0.1376731
## 3 121_TEMPEST_SW_H6_202311 22.0449 0.5923 142.0156 5.0936 0.6876138
## 4 131_TEMPEST_FW_E3_20231208 7.6841 0.2065 0.8009 0.0287 0.2396787
## 5 141_TEMPEST_C_C6_20231208 6.1506 0.1653 0.8842 0.0317 0.1918465
## 6 91_TEMPEST_FW_B4_20230807 4.9651 0.1334 32.8252 1.1773 0.1548690
##           Cl_mM      salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 0.027410437 0.0017818619 2.846054 7.797879e-05      50 0.001483
## 2 0.009726375 0.0006490536 2.805427 7.797879e-05      50 0.001483
## 3 4.006081805 0.2566481892 3.320780 7.797879e-05      50 0.001483
## 4 0.022592384 0.0014732263 2.834897 7.797879e-05      50 0.001483
## 5 0.024942172 0.0016237494 2.771666 7.797879e-05      50 0.001483
## 6 0.925957687 0.0593411364 2.919099 7.797879e-05      50 0.001483
## S04_Total_unspkd S04_Total_spkd S04_expctd_spkd spk_recovery S04_spks_flag
## 1 4.309859e-06 8.498318e-05 8.228865e-05 103.27449      YES
## 2 4.083385e-06 8.377006e-05 8.206217e-05 102.08121      YES
## 3 2.039463e-05 9.915848e-05 9.837342e-05 100.79805      YES
## 4 7.108871e-06 8.465003e-05 8.508766e-05 99.48567      YES
## 5 5.690168e-06 8.276194e-05 8.366896e-05 98.91594      YES
## 6 4.593414e-06 8.716428e-05 8.257220e-05 105.56129      YES
```

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

```
spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 YES           6      6     100
```

Pull out spikes and check with area calc

```
##Merge sample data with dilution data
Dilutions <- read.csv("Raw Data/TMP 2023 Sample Dilutions from Dionex.csv")
sampledat_dils <- merge(sampledat, Dilutions, by = "Sample_ID")
unique(sampledat_dils$Dilution)
```

```
## [1] 50.00 49.27
```



```
#Show me the data that we have from the calculations
head(sampledat_dils)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1 100_TEMPEST_C_B4_20230810 3.8752 0.1041 0.0705 0.0025 0.1208734
## 2 100_TEMPEST_C_B4_20230810_dup 4.2038 0.1129 0.0396 0.0014 0.1311229
## 3 101_TEMPEST_C_H6_20230807 4.6586 0.1252 0.9717 0.0349 0.1453088
## 4 101_TEMPEST_C_H6_20230807_spk 91.2445 2.4515 1.0566 0.0379 2.8460543
## 5 102_TEMPEST_C_H3_20230807 5.4909 0.1475 0.0087 0.0003 0.1712695
## 6 103_TEMPEST_C_F4_20230810 4.8364 0.1299 0.0371 0.0013 0.1508546
##           Cl_mM      salinity Dilution
## 1 0.0019887165 0.0001533935      50
## 2 0.0011170663 0.0000975572      50
## 3 0.0274104372 0.0017818619      50
## 4 0.0298053597 0.0019352762      50
## 5 0.0002454161 0.0000417209      50
## 6 0.0010465444 0.0000930397      50
```

```
#now use area of all samples and calculate undilution corrected concentrations in ug/mL
sampledat_dils$S04_ugmL <- ((sampledat_dils$S04_Area)-S04_Int)/S04_Slope
sampledat_dils$Cl_ugmL <- (sampledat_dils$Cl_Area-Cl_Int)/Cl_Slope
head(sampledat_dils)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1 100_TEMPEST_C_B4_20230810 3.8752 0.1041 0.0705 0.0025 0.1208734
## 2 100_TEMPEST_C_B4_20230810_dup 4.2038 0.1129 0.0396 0.0014 0.1311229
## 3 101_TEMPEST_C_H6_20230807 4.6586 0.1252 0.9717 0.0349 0.1453088
## 4 101_TEMPEST_C_H6_20230807_spk 91.2445 2.4515 1.0566 0.0379 2.8460543
## 5 102_TEMPEST_C_H3_20230807 5.4909 0.1475 0.0087 0.0003 0.1712695
## 6 103_TEMPEST_C_F4_20230810 4.8364 0.1299 0.0371 0.0013 0.1508546
##           Cl_mM      salinity Dilution   S04_ugmL   Cl_ugmL
## 1 0.0019887165 0.0001533935      50 0.07747505 0.0013990030
## 2 0.0011170663 0.0000975572      50 0.08402573 0.0007856160
## 3 0.0274104372 0.0017818619      50 0.09318179 0.0194660385
## 4 0.0298053597 0.0019352762      50 1.82486858 0.0211389122
## 5 0.0002454161 0.0000417209      50 0.10978180 0.0001722290
## 6 0.0010465444 0.0000930397      50 0.09668045 0.0007298536
```

```
#pull out any rows that have "spk" in the SampleID column
spks <- sampledat_dils %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1 101_TEMPEST_C_H6_20230807_spk 91.2445 2.4515 1.0566 0.0379 2.846054
## 2 111_TEMPEST_FW_I5_20231002_spk 89.9420 2.4165 0.3888 0.0139 2.805427
## 3 121_TEMPEST_SW_H6_202311_spk 106.4642 2.8604 141.0457 5.0588 3.320780
## 4 131_TEMPEST_FW_E3_20231208_spk 90.8868 2.4419 0.8722 0.0313 2.834897
## 5 141_TEMPEST_C_C6_20231208_spk 88.8596 2.3874 1.1352 0.0407 2.771666
## 6 91_TEMPEST_FW_B4_20230807_spk 93.5863 2.5144 0.0059 0.0002 2.919099
##           Cl_mM      salinity Dilution S04_ugmL   Cl_ugmL
## 1 0.0298053597 0.0019352762      50 1.824869 0.0211389122
```

```
## 2 0.0109675599 0.0007285616      50 1.798815 0.0077559229
## 3 3.9787221439 0.2548955799      50 2.129252 2.8209160213
## 4 0.0246036671 0.0016020654      50 1.817722 0.0174585901
## 5 0.0320225670 0.0020773064      50 1.777153 0.0227002609
## 6 0.0001664316 0.0000366613      50 1.871691 0.0001164666
```

```
## Only keep the columns that we need
spks <- spks[,c(1,9,10)]
head(spks)
```

```
##              Sample_ID Dilution S04_ugmL
## 1 101_TEMPEST_C_H6_20230807_spk      50 1.824869
## 2 111_TEMPEST_FW_I5_20231002_spk      50 1.798815
## 3 121_TEMPEST_SW_H6_202311_spk      50 2.129252
## 4 131_TEMPEST_FW_E3_20231208_spk      50 1.817722
## 5 141_TEMPEST_C_C6_20231208_spk      50 1.777153
## 6 91_TEMPEST_FW_B4_20230807_spk      50 1.871691
```

```
#remove the spk from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
colnames(spks) <- c('Sample_ID', 'Dilution', 'S04_ugmL_spk')
head(spks)
```

```
##              Sample_ID Dilution S04_ugmL_spk
## 1 101_TEMPEST_C_H6_20230807      50      1.824869
## 2 111_TEMPEST_FW_I5_20231002      50      1.798815
## 3 121_TEMPEST_SW_H6_202311      50      2.129252
## 4 131_TEMPEST_FW_E3_20231208      50      1.817722
## 5 141_TEMPEST_C_C6_20231208      50      1.777153
## 6 91_TEMPEST_FW_B4_20230807      50      1.871691
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat_dils, spks)
head(QAspks)
```

```
##              Sample_ID Dilution S04_ppm S04_Area  Cl_ppm Cl_Area
## 1 101_TEMPEST_C_H6_20230807      50 4.6586 0.1252 0.9717 0.0349
## 2 111_TEMPEST_FW_I5_20231002      50 4.4138 0.1186 0.3448 0.0124
## 3 121_TEMPEST_SW_H6_202311      50 22.0449 0.5923 142.0156 5.0936
## 4 131_TEMPEST_FW_E3_20231208      50 7.6841 0.2065 0.8009 0.0287
## 5 141_TEMPEST_C_C6_20231208      50 6.1506 0.1653 0.8842 0.0317
## 6 91_TEMPEST_FW_B4_20230807      50 4.9651 0.1334 32.8252 1.1773
##      S04_mM      Cl_mM      salinity      S04_ugmL      Cl_ugmL S04_ugmL_spk
## 1 0.1453088 0.027410437 0.0017818619 0.09318179 0.019466039      1.824869
## 2 0.1376731 0.009726375 0.0006490536 0.08826878 0.006919486      1.798815
## 3 0.6876138 4.006081805 0.2566481892 0.44088883 2.840321356      2.129252
## 4 0.2396787 0.022592384 0.0014732263 0.15370113 0.016008766      1.817722
## 5 0.1918465 0.024942172 0.0016237494 0.12303204 0.017681640      1.777153
## 6 0.1548690 0.925957687 0.0593411364 0.09928583 0.656496326      1.871691
```

```

#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250)      # in ug
spkvol <- 10          # in uL
spkvol <- spkvol/1000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol      # mmoles of SO4
head(QAspks)

```

```

##           Sample_ID Dilution SO4_ppm SO4_Area   Cl_ppm Cl_Area
## 1 101_TEMPEST_C_H6_20230807      50  4.6586   0.1252   0.9717  0.0349
## 2 111_TEMPEST_FW_I5_20231002      50  4.4138   0.1186   0.3448  0.0124
## 3 121_TEMPEST_SW_H6_202311      50 22.0449   0.5923 142.0156  5.0936
## 4 131_TEMPEST_FW_E3_20231208      50  7.6841   0.2065   0.8009  0.0287
## 5 141_TEMPEST_C_C6_20231208      50  6.1506   0.1653   0.8842  0.0317
## 6  91_TEMPEST_FW_B4_20230807      50  4.9651   0.1334  32.8252  1.1773
##      SO4_mM      Cl_mM      salinity   SO4_ugmL      Cl_ugmL SO4_ugmL_spk
## 1 0.1453088 0.027410437 0.0017818619 0.09318179 0.019466039      1.824869
## 2 0.1376731 0.009726375 0.0006490536 0.08826878 0.006919486      1.798815
## 3 0.6876138 4.006081805 0.2566481892 0.44088883 2.840321356      2.129252
## 4 0.2396787 0.022592384 0.0014732263 0.15370113 0.016008766      1.817722
## 5 0.1918465 0.024942172 0.0016237494 0.12303204 0.017681640      1.777153
## 6 0.1548690 0.925957687 0.0593411364 0.09928583 0.656496326      1.871691
##      SO4_spk_Conc
## 1              2.5
## 2              2.5
## 3              2.5
## 4              2.5
## 5              2.5
## 6              2.5

```

```

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(QAspks$Dilution == 50, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.85, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 51.16, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 49.27, 1483, QAspks$SampleVol)

```

```
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)
```

```
#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000
head(QAspks)
```

```
##           Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area
## 1 101_TEMPEST_C_H6_20230807      50  4.6586  0.1252  0.9717  0.0349
## 2 111_TEMPEST_FW_I5_20231002      50  4.4138  0.1186  0.3448  0.0124
## 3 121_TEMPEST_SW_H6_202311      50 22.0449  0.5923 142.0156  5.0936
## 4 131_TEMPEST_FW_E3_20231208      50  7.6841  0.2065  0.8009  0.0287
## 5 141_TEMPEST_C_C6_20231208      50  6.1506  0.1653  0.8842  0.0317
## 6 91_TEMPEST_FW_B4_20230807      50  4.9651  0.1334  32.8252  1.1773
##      S04_mM      Cl_mM      salinity      S04_ugmL      Cl_ugmL S04_ugmL_spk
## 1 0.1453088 0.027410437 0.0017818619 0.09318179 0.019466039      1.824869
## 2 0.1376731 0.009726375 0.0006490536 0.08826878 0.006919486      1.798815
## 3 0.6876138 4.006081805 0.2566481892 0.44088883 2.840321356      2.129252
## 4 0.2396787 0.022592384 0.0014732263 0.15370113 0.016008766      1.817722
## 5 0.1918465 0.024942172 0.0016237494 0.12303204 0.017681640      1.777153
## 6 0.1548690 0.925957687 0.0593411364 0.09928583 0.656496326      1.871691
##      S04_spk_Conc SampleVol
## 1          2.5          1.483
## 2          2.5          1.483
## 3          2.5          1.483
## 4          2.5          1.483
## 5          2.5          1.483
## 6          2.5          1.483
```

```
#gives us the total S04 in the sample in mmoles
```

```
QAspks$S04_Total_unspkd <- QAspks$S04_ugmL*QAspks$SampleVol
```

```
##total S04 in spiked sample in mmoles
```

```
QAspks$S04_Total_spkd <- (QAspks$S04_ugmL_spk)*(QAspks$SampleVol+spkvol)
```

```
QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
```

```
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')
```


```
head(QAspks)
```

```
##           Sample_ID Dilution S04_ppm S04_Area  Cl_ppm Cl_Area
## 1  101_TEMPEST_C_H6_20230807      50  4.6586   0.1252   0.9717  0.0349
## 2  111_TEMPEST_FW_I5_20231002      50  4.4138   0.1186   0.3448  0.0124
## 3   121_TEMPEST_SW_H6_202311      50 22.0449   0.5923 142.0156  5.0936
## 4  131_TEMPEST_FW_E3_20231208      50  7.6841   0.2065   0.8009  0.0287
## 5  141_TEMPEST_C_C6_20231208      50  6.1506   0.1653   0.8842  0.0317
## 6   91_TEMPEST_FW_B4_20230807      50  4.9651   0.1334  32.8252  1.1773
##      S04_mM      Cl_mM      salinity  S04_ugmL      Cl_ugmL S04_ugmL_spk
## 1 0.1453088 0.027410437 0.0017818619 0.09318179 0.019466039      1.824869
## 2 0.1376731 0.009726375 0.0006490536 0.08826878 0.006919486      1.798815
## 3 0.6876138 4.006081805 0.2566481892 0.44088883 2.840321356      2.129252
## 4 0.2396787 0.022592384 0.0014732263 0.15370113 0.016008766      1.817722
## 5 0.1918465 0.024942172 0.0016237494 0.12303204 0.017681640      1.777153
## 6 0.1548690 0.925957687 0.0593411364 0.09928583 0.656496326      1.871691
##      S04_spk_Conc SampleVol S04_Total_unspkd S04_Total_spkd S04_expctd_spkd
## 1           2.5         1.483         0.1381886         2.724529         2.638189
## 2           2.5         1.483         0.1309026         2.685630         2.630903
## 3           2.5         1.483         0.6538381         3.178973         3.153838
## 4           2.5         1.483         0.2279388         2.713860         2.727939
## 5           2.5         1.483         0.1824565         2.653289         2.682457
## 6           2.5         1.483         0.1472409         2.794435         2.647241
##      spk_recovery S04_spks_flag
## 1      103.27271             YES
## 2      102.08019             YES
## 3      100.79696             YES
## 4       99.48389             YES
## 5       98.91266             YES
## 6      105.56027             YES
```

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

```
spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 YES             6     6    100
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_ ', fixed=TRUE)))
colnames(IDs) <- c("Number", "Project", "Plot", "Grid", "Date")
head(IDs)
```

```
##   Number Project Plot Grid   Date
## 1    100 TEMPEST   C   B4 20230810
```

```
## 2    101 TEMPEST    C    H6 20230807
## 3    102 TEMPEST    C    H3 20230807
## 4    103 TEMPEST    C    F4 20230810
## 5    104 TEMPEST    C    C6 20230807
## 6    105 TEMPEST   SW    I5 20231002
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)

alldat$Number <- NULL

alldat$Sample_ID <- paste(alldat$Project, alldat$Plot, alldat$Grid, alldat$Date, sep = "_")

head(alldat)
```

```
##   Project Plot Grid    Date      Sample_ID  SO4_ppm SO4_Area  Cl_ppm
## 1 TEMPEST    C    B4 20230810 TEMPEST_C_B4_20230810  3.8752  0.1041  0.0705
## 2 TEMPEST    C    H6 20230807 TEMPEST_C_H6_20230807  4.6586  0.1252  0.9717
## 3 TEMPEST    C    H3 20230807 TEMPEST_C_H3_20230807  5.4909  0.1475  0.0087
## 4 TEMPEST    C    F4 20230810 TEMPEST_C_F4_20230810  4.8364  0.1299  0.0371
## 5 TEMPEST    C    C6 20230807 TEMPEST_C_C6_20230807  4.8123  0.1293  3.7102
## 6 TEMPEST   SW    I5 20231002 TEMPEST_SW_I5_20231002 114.6553  3.0805 281.4327
##   Cl_Area  SO4_mM      Cl_mM    salinity
## 1  0.0025 0.1208734 0.0019887165 0.0001533935
## 2  0.0349 0.1453088 0.0274104372 0.0017818619
## 3  0.0003 0.1712695 0.0002454161 0.0000417209
## 4  0.0013 0.1508546 0.0010465444 0.0000930397
## 5  0.1331 0.1501029 0.1046600846 0.0067303314
## 6 10.0940 3.5762726 7.9388631876 0.5085748889
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_TEMPEST_Dionex_ClSo4_Final_2023a.csv")
```

#Change

END