

Dionex_COMPASS_TEMPEST

Stephanie J. Wilson

2023-06-06

Daily Set up

Read in and Format the raw data - change wd & file names (Make sure files are UTF-8 encoded)

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_TEMPEST_20230606_S04a.txt", sep='\t' , header=T, skip=3)
head(Sdat)
```

```
##      X              X.1      X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 44 TEMPEST_C_C6_20230603 Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 2 45 TEMPEST_C_H6_20230603 Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 3 46 TEMPEST_C_I5_20230603 Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
## 4 47 TEMPEST_C_H6_20230605 Unknown  4.373    2.3894    0.47    0.0544    0.35
## 5 48 TEMPEST_FW_C3_20230603 Unknown  4.383    6.9627    0.69    0.1403    0.90
## 6 49 TEMPEST_FW_D5_20230603 Unknown  n.a.      n.a.      n.a.      n.a.      n.a.
##      IC.S04.5
## 1      n.a.
## 2      n.a.
## 3      n.a.
## 4      M
## 5      M
## 6      n.a.
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[ ,c(2,5,7)] # dont need this here
head(Sdat)
```

```
##              X.1 IC.S04.1 IC.S04.3
## 1 TEMPEST_C_C6_20230603  n.a.      n.a.
## 2 TEMPEST_C_H6_20230603  n.a.      n.a.
## 3 TEMPEST_C_I5_20230603  n.a.      n.a.
## 4 TEMPEST_C_H6_20230605  2.3894    0.0544
## 5 TEMPEST_FW_C3_20230603  6.9627    0.1403
## 6 TEMPEST_FW_D5_20230603  n.a.      n.a.
```

```
## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "S04_ppm", "S04_Area")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat$S04_Area <- as.numeric(Sdat$S04_Area)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##           Sample_ID S04_ppm S04_Area
## 1 TEMPEST_C_C6_20230603      NA      NA
## 2 TEMPEST_C_H6_20230603      NA      NA
## 3 TEMPEST_C_I5_20230603      NA      NA
## 4 TEMPEST_C_H6_20230605  2.3894  0.0544
## 5 TEMPEST_FW_C3_20230603  6.9627  0.1403
## 6 TEMPEST_FW_D5_20230603      NA      NA
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("Raw Data/COMPASS_TEMPEST_20230606_Cla.txt",sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X              X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 44 TEMPEST_C_C6_20230603 Unknown n.a.    n.a.    n.a.    n.a.    n.a.
## 2 45 TEMPEST_C_H6_20230603 Unknown n.a.    n.a.    n.a.    n.a.    n.a.
## 3 46 TEMPEST_C_I5_20230603 Unknown n.a.    n.a.    n.a.    n.a.    n.a.
## 4 47 TEMPEST_C_H6_20230605 Unknown n.a.    n.a.    n.a.    n.a.    n.a.
## 5 48 TEMPEST_FW_C3_20230603 Unknown n.a.    n.a.    n.a.    n.a.    n.a.
## 6 49 TEMPEST_FW_D5_20230603 Unknown n.a.    n.a.    n.a.    n.a.    n.a.
##      IC.Cl.5
## 1      n.a.
## 2      n.a.
## 3      n.a.
## 4      n.a.
## 5      n.a.
## 6      n.a.
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[ ,c(2,5,7)]
head(Cldat)
```

```
##           X.1 IC.Cl.1 IC.Cl.3
## 1 TEMPEST_C_C6_20230603  n.a.  n.a.
## 2 TEMPEST_C_H6_20230603  n.a.  n.a.
## 3 TEMPEST_C_I5_20230603  n.a.  n.a.
## 4 TEMPEST_C_H6_20230605  n.a.  n.a.
## 5 TEMPEST_FW_C3_20230603  n.a.  n.a.
## 6 TEMPEST_FW_D5_20230603  n.a.  n.a.
```

```
## Name the columns correctly
colnames(Cldat) <- c( "Sample_ID", "Cl_ppm", "Cl_Area")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat$Cl_Area <- as.numeric(Cldat$Cl_Area)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##           Sample_ID Cl_ppm Cl_Area
## 1 TEMPEST_C_C6_20230603    NA     NA
## 2 TEMPEST_C_H6_20230603    NA     NA
## 3 TEMPEST_C_I5_20230603    NA     NA
## 4 TEMPEST_C_H6_20230605    NA     NA
## 5 TEMPEST_FW_C3_20230603    NA     NA
## 6 TEMPEST_FW_D5_20230603    NA     NA
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##   Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1          168.9031  23.3576 203.7164 279.6752
## 2          168.9031  23.3576  29.8315  32.6201
## 3          168.9031  23.3576   0.0014   0.0001
## 4          168.9031  23.3576  40.9777  53.7128
## 5          168.9031  23.3576        NA        NA
## 6           13.6022   2.4364 203.7164 279.6752
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##   Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 26 Lab Blank    NA        NA    NA    NA
## 27 Lab Blank    NA        NA 0.0079 0.0108
## 28 Lab Blank    NA        NA    NA    NA
## 29 Lab Blank    NA        NA    NA    NA
## 30 Lab Blank    NA        NA 0.0114 0.0156
## 31 Lab Blank    NA        NA    NA    NA
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 107 Standard 1  0.4622  0.5264  5.2376  7.1905
## 108 Standard 1  0.4622  0.5264  5.1529  7.0742
## 109 Standard 1  0.4374  0.4982  5.2376  7.1905
## 110 Standard 1  0.4374  0.4982  5.1529  7.0742
## 111 Standard 2  1.0048  1.1446 11.4708 15.7479
## 112 Standard 2  1.0048  1.1446 11.3597 15.5954
```

```
all_dat[is.na(all_dat)] <- 0
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grepl("Standard", all_dat$Sample_ID),]
head(stds)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 107 Standard 1  0.4622  0.5264  5.2376  7.1905
## 108 Standard 1  0.4622  0.5264  5.1529  7.0742
## 109 Standard 1  0.4374  0.4982  5.2376  7.1905
## 110 Standard 1  0.4374  0.4982  5.1529  7.0742
## 111 Standard 2  1.0048  1.1446 11.4708 15.7479
## 112 Standard 2  1.0048  1.1446 11.3597 15.5954
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(S04_ppm), sd = sd(S04_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean      sd      cv flag
##   <fct>      <dbl>   <dbl> <dbl> <chr>
## 1 Standard 1  0.450  0.0143  3.18 NO, rerun
## 2 Standard 2  0.993  0.0138  1.39 YES
## 3 Standard 3  2.02   0.0309  1.53 YES
## 4 Standard 4 10.5    NA      NA    <NA>
## 5 Standard 5 20.5    NA      NA    <NA>
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean      sd      cv flag
##   <fct>      <dbl>   <dbl> <dbl> <chr>
## 1 Standard 1  5.20  0.0489  0.941 YES
```

```
## 2 Standard 2 11.4 0.0641 0.562 YES
## 3 Standard 3 20.8 0.117 0.564 YES
## 4 Standard 4 103. NA NA <NA>
## 5 Standard 5 204. NA NA <NA>
```

```
lmS <- lm(stds$S04_Area ~ stds$S04_ppm)
S04_sum <- summary(lmS)
S04_Slope <- S04_sum$coefficients[2, 1]
S04_Int <- S04_sum$coefficients[1, 1]
```

```
lmCl <- lm(stds$Cl_Area ~ stds$Cl_ppm)
Cl_sum <- summary(lmCl)
Cl_Slope <- Cl_sum$coefficients[2, 1]
Cl_Int <- Cl_sum$coefficients[1, 1]
```

Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 126 TEMPEST_C_C6_20230515 0.3224 0.3338 30.5616 38.1427
## 127 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000
## 128 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000
## 129 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000
## 130 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000
## 131 TEMPEST_C_F6_20230515 0.0000 0.0000 28.2531 35.2615
```

```
# Constants needed for calculations:
clmw <- 35.45 #molecular weight of Chloride: 35.45
smw <- 32.06 #molecular weight of sulfur: 32.06
```

```
# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
# = ((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
```

```
## 126 TEMPEST_C_C6_20230515 0.3224 0.3338 30.5616 38.1427 0.01005614 0.8621044
## 127 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 128 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 129 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 130 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 131 TEMPEST_C_F6_20230515 0.0000 0.0000 28.2531 35.2615 0.00000000 0.7969845
##      salinity
## 126 0.05525081
## 127 0.00002600
## 128 0.00002600
## 129 0.00002600
## 130 0.00002600
## 131 0.05107935
```

Pull out dups and check with percent difference

```
#Show me the data that we have from the calculations
head(sampledat)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 126 TEMPEST_C_C6_20230515 0.3224 0.3338 30.5616 38.1427 0.01005614 0.8621044
## 127 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 128 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 129 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 130 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 131 TEMPEST_C_F6_20230515 0.0000 0.0000 28.2531 35.2615 0.00000000 0.7969845
##      salinity
## 126 0.05525081
## 127 0.00002600
## 128 0.00002600
## 129 0.00002600
## 130 0.00002600
## 131 0.05107935
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "_dup"))      #have to change this to match data
head(dups)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM
## 1 TEMPEST_C_I5_20230518_dup 4.1578 4.3055 25.5009 31.8266 0.12968808
## 2 TEMPEST_SW_H3_20230603_dup 0.0000 0.0000 0.0000 0.0000 0.00000000
## 3 TEMPEST_SW_H6_20230515_dup 0.4053 0.4197 36.6535 45.7458 0.01264192
##      Cl_mM      salinity
## 1 0.7193484 0.04610613
## 2 0.0000000 0.00002600
## 3 1.0339492 0.06625887
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "_dup")) %>%
```

```
filter(!str_detect(Sample_ID, "_spk"))
head(sampled2)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_C6_20230515  0.3224  0.3338 30.5616 38.1427 0.01005614 0.8621044
## 2 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 3 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 4 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 5 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 6 TEMPEST_C_F6_20230515  0.0000  0.0000 28.2531 35.2615 0.00000000 0.7969845
##      salinity
## 1 0.05525081
## 2 0.00002600
## 3 0.00002600
## 4 0.00002600
## 5 0.00002600
## 6 0.05107935
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[ , -c(2:5)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', "salinity_dup")
head(dups)
```

```
##           Sample_ID S04_mM_dup Cl_mM_dup salinity_dup
## 1 TEMPEST_C_I5_20230518 0.12968808 0.7193484  0.04610613
## 2 TEMPEST_SW_H3_20230603 0.00000000 0.0000000  0.00002600
## 3 TEMPEST_SW_H6_20230515 0.01264192 1.0339492  0.06625887
```

```
#put it back together with the old data set and look for duplicates
QAdups <- merge(sampled2, dups)
head(QAdups)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_I5_20230518  4.4207  4.5778 25.4606 31.7763 0.13788833 0.7182116
## 2 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 3 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 4 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 5 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 6 TEMPEST_SW_H6_20230515  0.3839  0.3975 28.1246 35.1012 0.01197442 0.7933597
##      salinity S04_mM_dup Cl_mM_dup salinity_dup
## 1 0.04603330 0.12968808 0.7193484  0.04610613
## 2 0.00002600 0.00000000 0.0000000  0.00002600
## 3 0.00002600 0.00000000 0.0000000  0.00002600
## 4 0.00002600 0.00000000 0.0000000  0.00002600
## 5 0.00002600 0.00000000 0.0000000  0.00002600
## 6 0.05084715 0.01264192 1.0339492  0.06625887
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*10
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')

head(QAdups)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_I5_20230518  4.4207  4.5778 25.4606 31.7763 0.13788833 0.7182116
## 2 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 3 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 4 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 5 TEMPEST_SW_H3_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 6 TEMPEST_SW_H6_20230515  0.3839  0.3975 28.1246 35.1012 0.01197442 0.7933597
##      salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 0.04603330 0.12968808 0.7193484  0.04610613      6.129277      YES
## 2 0.00002600 0.00000000 0.0000000  0.00002600      NaN      <NA>
## 3 0.00002600 0.00000000 0.0000000  0.00002600      NaN      <NA>
## 4 0.00002600 0.00000000 0.0000000  0.00002600      NaN      <NA>
## 5 0.00002600 0.00000000 0.0000000  0.00002600      NaN      <NA>
## 6 0.05084715 0.01264192 1.0339492  0.06625887      5.423213      YES
##      Cl_dups_chk Cl_dups_flag
## 1  0.1581586      YES
## 2      NaN      <NA>
## 3      NaN      <NA>
## 4      NaN      <NA>
## 5      NaN      <NA>
## 6 26.3326649      NO, rerun
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```


```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```

```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_bar()').
```



```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_bar()').
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups <- Perc_dups %>%
  filter(!is.na(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
Perc_dups1 <- Perc_dups1 %>%
  filter(!is.na(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows      Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         2 NO, rerun          1     6  33.33333  16.66667
## 2  YES         2   YES              1     6  33.33333  16.66667
```

Pull out dups and check with cv

```
#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)
```


```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_I5_20230518 4.4207 4.5778 25.4606 31.7763 0.13788833 0.7182116
## 2 TEMPEST_SW_H3_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 3 TEMPEST_SW_H3_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 4 TEMPEST_SW_H3_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 5 TEMPEST_SW_H3_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 6 TEMPEST_SW_H6_20230515 0.3839 0.3975 28.1246 35.1012 0.01197442 0.7933597
##      salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk S04_dups_flag
## 1 0.04603330 0.12968808 0.7193484 0.04610613      6.129277      YES
## 2 0.00002600 0.00000000 0.0000000 0.00002600      NaN      <NA>
## 3 0.00002600 0.00000000 0.0000000 0.00002600      NaN      <NA>
## 4 0.00002600 0.00000000 0.0000000 0.00002600      NaN      <NA>
## 5 0.00002600 0.00000000 0.0000000 0.00002600      NaN      <NA>
## 6 0.05084715 0.01264192 1.0339492 0.06625887      5.423213      YES
##      Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1 0.1581586      YES      4.334053      YES
## 2      NaN      <NA>      NaN      <NA>
## 3      NaN      <NA>      NaN      <NA>
## 4      NaN      <NA>      NaN      <NA>
## 5      NaN      <NA>      NaN      <NA>
## 6 26.3326649      NO, rerun      3.834791      YES
```

```
#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)
```

```
ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```

```
## Warning: Removed 4 rows containing missing values or values outside the scale range
## ('geom_bar()').
## Removed 4 rows containing missing values or values outside the scale range
## ('geom_bar()').
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups <- Perc_dups %>%
  filter(!is.na(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
Perc_dups1 <- Perc_dups1 %>%
  filter(!is.na(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
```

```
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows      Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         2 NO, rerun          1     6  33.33333  16.66667
## 2  YES         2    YES          1     6  33.33333  16.66667
```

Pull out spikes and check with dionex output conc.

```
##Merge sample data with dilution data
Dilutions <- read.csv("Raw Data/TMP 2023 Sample Dilutions from Dionex.csv")
sampledat_dils <- merge(sampledat, Dilutions, by = "Sample_ID")
unique(sampledat_dils$Dilution)
```

```
## [1] 1.10 50.00 5.15 56.51 4.53
```

```
#Show me the data that we have from the calculations
head(sampledat_dils)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_C6_20230515 0.3224 0.3338 30.5616 38.1427 0.01005614 0.8621044
## 2 TEMPEST_C_C6_20230515 0.3224 0.3338 30.5616 38.1427 0.01005614 0.8621044
## 3 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 4 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 5 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
## 6 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000 0.0000 0.00000000 0.0000000
##      salinity Dilution
## 1 0.05525081      1.10
## 2 0.05525081     50.00
## 3 0.00002600      5.15
## 4 0.00002600     56.51
## 5 0.00002600      5.15
## 6 0.00002600     56.51
```

```
#pull out any rows that have "spk" in the SampleID column
spks <- sampledat_dils %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM Cl_mM
## 1 TEMPEST_SW_H3_20230603_spk 168.9031 3.4046      0      0 5.268344      0
##      salinity Dilution
## 1 2.6e-05     56.51
```

```
#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[,c(1,6)]
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
head(spks)
```

```
##           Sample_ID S04_mM_spk
## 1 TEMPEST_SW_H3_20230603 5.268344
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledats, spks)
head(QAspks)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM salinity
## 1 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 2 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 3 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 4 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 5 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 6 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
##   Dilution S04_mM_spk
## 1    56.51 5.268344
## 2     5.15 5.268344
## 3    56.51 5.268344
## 4    56.51 5.268344
## 5     5.15 5.268344
## 6    56.51 5.268344
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10              # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol      # mmoles of S04
head(QAspks)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM salinity
## 1 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 2 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 3 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 4 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 5 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 6 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
##   Dilution S04_mM_spk S04_spk_Conc
## 1    56.51 5.268344 7.797879e-05
## 2     5.15 5.268344 7.797879e-05
## 3    56.51 5.268344 7.797879e-05
## 4    56.51 5.268344 7.797879e-05
## 5     5.15 5.268344 7.797879e-05
## 6    56.51 5.268344 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
# QAspks$Dilution <- 1
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
```

```

# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
unique(QAspks$Dilution)

```

```
## [1] 56.51 5.15
```

```

QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(QAspks$Dilution == 50, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.85, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 51.16, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 49.27, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.51, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 5.15, 1500, QAspks$SampleVol)

# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 1483, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

```

```

##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM salinity
## 1 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 2 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 3 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 4 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 5 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 6 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
##   Dilution S04_mM_spk S04_spk_Conc SampleVol
## 1    56.51   5.268344 7.797879e-05   0.001483
## 2     5.15   5.268344 7.797879e-05   0.001500

```

```
## 3      56.51      5.268344 7.797879e-05 0.001483
## 4      56.51      5.268344 7.797879e-05 0.001483
## 5       5.15      5.268344 7.797879e-05 0.001500
## 6      56.51      5.268344 7.797879e-05 0.001483
```

```
#gives us the total SO4 in the sample in mmoles
```

```
QAspks$SO4_Total_unspkd <- (QAspks$SO4_mM/QAspks$Dilution)*(QAspks$SampleVol)
```

```
##total SO4 in spiked sample in mmoles
```

```
QAspks$SO4_Total_spkd <- (QAspks$SO4_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)
```

```
QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$SO4_Total_spkd/QAspks$SO4_expctd_spkd)*100
```

```
QAspks$SO4_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun
```


```
head(QAspks)
```

```
##           Sample_ID SO4_ppm SO4_Area Cl_ppm Cl_Area SO4_mM Cl_mM salinity
## 1 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 2 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 3 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 4 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 5 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
## 6 TEMPEST_SW_H3_20230603      0      0      0      0      0      0 2.6e-05
##   Dilution SO4_mM_spk SO4_spk_Conc SampleVol SO4_Total_unspkd SO4_Total_spkd
## 1    56.51    5.268344 7.797879e-05 0.001483              0 0.0001391902
## 2     5.15    5.268344 7.797879e-05 0.001500              0 0.0015446988
## 3    56.51    5.268344 7.797879e-05 0.001483              0 0.0001391902
## 4    56.51    5.268344 7.797879e-05 0.001483              0 0.0001391902
## 5     5.15    5.268344 7.797879e-05 0.001500              0 0.0015446988
## 6    56.51    5.268344 7.797879e-05 0.001483              0 0.0001391902
##   SO4_expctd_spkd spk_recovery SO4_spks_flag
## 1 7.797879e-05      178.4975      NO, rerun
## 2 7.797879e-05     1980.9218      NO, rerun
## 3 7.797879e-05      178.4975      NO, rerun
## 4 7.797879e-05      178.4975      NO, rerun
## 5 7.797879e-05     1980.9218      NO, rerun
## 6 7.797879e-05      178.4975      NO, rerun
```

```
#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work o
```

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

```
spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun           8     8    100
```

Pull out spikes and check with area calc

```
##Merge sample data with dilution data
Dilutions <- read.csv("Raw Data/TMP 2023 Sample Dilutions from Dionex.csv")
sampledat_dils <- merge(sampledat, Dilutions, by = "Sample_ID")
unique(sampledat_dils$Dilution)
```

```
## [1] 1.10 50.00 5.15 56.51 4.53
```



```
#Show me the data that we have from the calculations
head(sampledats_dils)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_C6_20230515  0.3224  0.3338 30.5616 38.1427 0.01005614 0.8621044
## 2 TEMPEST_C_C6_20230515  0.3224  0.3338 30.5616 38.1427 0.01005614 0.8621044
## 3 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 4 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 5 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 6 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
##      salinity Dilution
## 1 0.05525081      1.10
## 2 0.05525081     50.00
## 3 0.00002600      5.15
## 4 0.00002600     56.51
## 5 0.00002600      5.15
## 6 0.00002600     56.51
```

```
#now use area of all samples and calculate undilution corrected concentrations in ug/mL
sampledats_dils$S04_ugmL <- ((sampledats_dils$S04_Area)-S04_Int)/S04_Slope
sampledats_dils$Cl_ugmL <- (sampledats_dils$Cl_Area-Cl_Int)/Cl_Slope
head(sampledats_dils)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_C_C6_20230515  0.3224  0.3338 30.5616 38.1427 0.01005614 0.8621044
## 2 TEMPEST_C_C6_20230515  0.3224  0.3338 30.5616 38.1427 0.01005614 0.8621044
## 3 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 4 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 5 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
## 6 TEMPEST_C_C6_20230603  0.0000  0.0000  0.0000  0.0000 0.00000000 0.0000000
##      salinity Dilution      S04_ugmL      Cl_ugmL
## 1 0.05525081      1.10 2.93060e-01 2.778329e+01
## 2 0.05525081     50.00 2.93060e-01 2.778329e+01
## 3 0.00002600      5.15 1.61704e-05 7.223654e-06
## 4 0.00002600     56.51 1.61704e-05 7.223654e-06
## 5 0.00002600      5.15 1.61704e-05 7.223654e-06
## 6 0.00002600     56.51 1.61704e-05 7.223654e-06
```

```
#pull out any rows that have "spk" in the SampleID column
spks <- sampledats_dils %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area      S04_mM      Cl_mM
## 1 TEMPEST_SW_H3_20230603_spk 168.9031  3.4046      0      0 5.268344      0
##      salinity Dilution S04_ugmL      Cl_ugmL
## 1 2.6e-05     56.51 2.988922 7.223654e-06
```

```
## Only keep the columns that we need
spks <- spks[,c(1,9,10)]
head(spks)
```

```
##           Sample_ID Dilution S04_ugmL
## 1 TEMPEST_SW_H3_20230603_spk      56.51 2.988922
```

```
#remove the spk from these IDs so we will have duplicate sample names
```

```
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
colnames(spks) <- c('Sample_ID', 'Dilution', 'S04_ugmL_spk')
head(spks)
```

```
##           Sample_ID Dilution S04_ugmL_spk
## 1 TEMPEST_SW_H3_20230603      56.51      2.988922
```

```
#put it back together with the old data set and look for duplicates
```

```
QAspks <- merge(sampledat_dils, spks)
head(QAspks)
```

```
##           Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM
## 1 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
## 2 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
## 3 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
## 4 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
##   salinity   S04_ugmL   Cl_ugmL S04_ugmL_spk
## 1  2.6e-05 1.61704e-05 7.223654e-06      2.988922
## 2  2.6e-05 1.61704e-05 7.223654e-06      2.988922
## 3  2.6e-05 1.61704e-05 7.223654e-06      2.988922
## 4  2.6e-05 1.61704e-05 7.223654e-06      2.988922
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
```

```
spkconcc <- (250)      # in ug
spkvol <- 10           # in uL
spkvol <- spkvol/1000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconcc)*spkvol      # mmoles of S04
head(QAspks)
```

```
##           Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM
## 1 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
## 2 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
## 3 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
## 4 TEMPEST_SW_H3_20230603      56.51      0      0      0      0      0      0
##   salinity   S04_ugmL   Cl_ugmL S04_ugmL_spk S04_spk_Conc
## 1  2.6e-05 1.61704e-05 7.223654e-06      2.988922      2.5
## 2  2.6e-05 1.61704e-05 7.223654e-06      2.988922      2.5
## 3  2.6e-05 1.61704e-05 7.223654e-06      2.988922      2.5
## 4  2.6e-05 1.61704e-05 7.223654e-06      2.988922      2.5
```

```
#need to determine dilution factors and initial amount of sample added
```

```
#if your samples are all the same dilution just use the first line
```

```
#for Steph / COMPASS this depends on the site so...
```

```
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
```

```

##Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
##Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(QAspks$Dilution == 50, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 1.1, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 5.15, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 4.53, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.51, 1483, QAspks$SampleVol)

# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000
head(QAspks)

```

```

##          Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM
## 1 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
## 2 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
## 3 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
## 4 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
##   salinity   S04_ugmL   Cl_ugmL S04_ugmL_spk S04_spk_Conc SampleVol
## 1  2.6e-05 1.61704e-05 7.223654e-06    2.988922        2.5    1.483
## 2  2.6e-05 1.61704e-05 7.223654e-06    2.988922        2.5    1.483
## 3  2.6e-05 1.61704e-05 7.223654e-06    2.988922        2.5    1.483
## 4  2.6e-05 1.61704e-05 7.223654e-06    2.988922        2.5    1.483

```

```

#gives us the total SO4 in the sample in mmoles
QAspks$SO4_Total_unspkd <- QAspks$S04_ugmL*QAspks$SampleVol

```

```

##total S04 in spiked sample in mmoles
QAspks$S04_Total_spkd <- (QAspks$S04_ugmL_spk)*(QAspks$SampleVol+spkvol)

QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')


head(QAspks)

##           Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area S04_mM Cl_mM
## 1 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
## 2 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
## 3 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
## 4 TEMPEST_SW_H3_20230603    56.51      0      0      0      0      0      0
##   salinity   S04_ugmL   Cl_ugmL S04_ugmL_spk S04_spk_Conc SampleVol
## 1  2.6e-05 1.61704e-05 7.223654e-06    2.988922         2.5     1.483
## 2  2.6e-05 1.61704e-05 7.223654e-06    2.988922         2.5     1.483
## 3  2.6e-05 1.61704e-05 7.223654e-06    2.988922         2.5     1.483
## 4  2.6e-05 1.61704e-05 7.223654e-06    2.988922         2.5     1.483
##   S04_Total_unspkd S04_Total_spkd S04_expctd_spkd spk_recovery S04_spks_flag
## 1    2.398071e-05    4.462461    2.500024    178.4967    NO, rerun
## 2    2.398071e-05    4.462461    2.500024    178.4967    NO, rerun
## 3    2.398071e-05    4.462461    2.500024    178.4967    NO, rerun
## 4    2.398071e-05    4.462461    2.500024    178.4967    NO, rerun

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)

spksbar

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 NO, rerun      4      4     100
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_ ', fixed=TRUE)))
colnames(IDs) <- c("Project" , "Plot", "Grid", "Date")
head(IDs)
```

```
##   Project Plot Grid   Date
## 1 TEMPEST   C   C6 20230515
```

```
## 2 TEMPEST C C6 20230603
## 3 TEMPEST C C6 20230603
## 4 TEMPEST C C6 20230603
## 5 TEMPEST C C6 20230603
## 6 TEMPEST C F6 20230515
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)
head(alldat)
```

```
## Project Plot Grid Date Sample_ID S04_ppm S04_Area Cl_ppm
## 1 TEMPEST C C6 20230515 TEMPEST_C_C6_20230515 0.3224 0.3338 30.5616
## 2 TEMPEST C C6 20230603 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000
## 3 TEMPEST C C6 20230603 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000
## 4 TEMPEST C C6 20230603 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000
## 5 TEMPEST C C6 20230603 TEMPEST_C_C6_20230603 0.0000 0.0000 0.0000
## 6 TEMPEST C F6 20230515 TEMPEST_C_F6_20230515 0.0000 0.0000 28.2531
## Cl_Area S04_mM Cl_mM salinity
## 1 38.1427 0.01005614 0.8621044 0.05525081
## 2 0.0000 0.00000000 0.0000000 0.00002600
## 3 0.0000 0.00000000 0.0000000 0.00002600
## 4 0.0000 0.00000000 0.0000000 0.00002600
## 5 0.0000 0.00000000 0.0000000 0.00002600
## 6 35.2615 0.00000000 0.7969845 0.05107935
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_TEMPEST_Dionex_ClSo4_Final_20230606.csv")
```

#Cha

END