

# COMPASS: TEMPEST Discrete DOC Data QAQC

September 2025

2025-10-02

## Run Information

```
#identify which section you are in  
cat("Run Information")
```

```
## Run Information
```

```
#a link to the Gitbook or whatever protocol you are using for this analysis  
#steph will add this soon
```

```
#anything that needs to be changed do this in the first chunk
```

```
Date_Run = "07/03/25"
```

```
Run_by = "Stephanie J. Wilson"
```

```
Script_run_by = "Stephanie J. Wilson"
```

```
run_notes = "Run by Watershed Science at SERC, no duplicates included"
```

```
#file path and name for summary file
```

```
raw_file_name = "tmp_doc_raw_data_2025/TMP_202509.txt"
```

```
#file path and name for the all peaks file
```

```
raw_allpeaks_name = "tmp_doc_raw_data_2025/TMP_202509_allpeaks.txt"
```

```
#file path and name for processed data after QAQC
```

```
processed_file_name = "tmp_doc_processed_data_2025/TMP_PW_DOC_Processed_202509.csv"
```

```
#check standard concentrations - Update if running different checks:
```

```
chk_std_c = 30
```

```
chk_std_n = 2.5
```

```
#Log path
```

```
Log_path = "tmp_doc_raw_data_2025/COMPASS_TMP_TOCTN_QAQClog_2025.csv"
```

## Setup

### Pull in active porewater tracking inventory sheet

```
## File already exists. No download needed.
```

## Import Data Functions

### Import Sample Data

```
## Import Sample Data
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## # A tibble: 3 x 4
```

```
##   sample_name      npoc_raw tdn_raw run_datetime  
##   <chr>          <dbl>   <dbl> <chr>  
## 1 TMP_SW_D5_20250926    15.7    0.605 10/1/2025 10:45:08 PM  
## 2 TMP_SW_H3_20250929    17.7    0.643 10/1/2025 11:05:35 PM  
## 3 TMP_SW_H6_20250926    20.8    0.868 10/1/2025 11:33:43 PM
```


## Assessing standard Curves

```
## Assess the Standard Curve
```


```
## New names:
```


```
## 'geom_smooth()' using formula = 'y ~ x'
```

```
## * '' -> '...18'
```



```
## 'geom_smooth()' using formula = 'y ~ x'
```





```
## [1] "NPOC Curve r2 GOOD"
```

```
## [1] "TN Curve r2 GOOD"
```

## Assess Check Standards


```
## Assess the Check Standards
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## [1] "Carbon Check Standard RSD within Range"
```

```
## [1] "Nitrogen Check Standard RSD within Range"
```



```
## [1] ">60% of Carbon Check Standards are within range of the expected concentration"
```

```
## [1] ">60% of Nitrogen Check Standards are within range of the expected concentration"
```

## Assess Blanks


```
## Assess Blanks
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## [1] ">60% of Carbon Blank concentrations are below the lower 25% quartile of samples"
```

```
## [1] ">60% of Nitrogen Blank concentrations are below the lower 25% quartile of samples"
```



ink Conc <25% Quartile Samples  NO, r Blank Conc <25% Quartile Samples  NO, rerun

```
## carbon blanks:
```


```
## [1] 0.7973995
```

```
## nitrogen blanks:
```


```
## [1] 0.03649595
```

## Sample Flagging

## Sample Flagging









Visualize Data by Plot

## Visualize Data

##	Site_Code	Plot	Grid_Square	Date
## 1	TMP	SW	D5	20250926
## 2	TMP	SW	H3	20250929
## 3	TMP	SW	H6	20250926

##	Site_Code	Plot	Grid_Square	Date	sample_name	npoc_raw	tdn_raw
## 1	TMP	SW	D5	20250926	TMP_SW_D5_20250926	15.73	0.6051
## 2	TMP	SW	H3	20250929	TMP_SW_H3_20250929	17.70	0.6428
## 3	TMP	SW	H6	20250926	TMP_SW_H6_20250926	20.77	0.8676
##	run_datetime		npoc_flag	tdn_flag			
## 1	10/1/2025	10:45:08 PM					
## 2	10/1/2025	11:05:35 PM					
## 3	10/1/2025	11:33:43 PM					





## Convert data from mg/L to uMoles/L

### Add in/check metadata

```
## Check Sample IDs with Metadata
```

```
## # A tibble: 3 x 2
##   sample_name      metadata_recorded
##   <chr>            <lgl>
## 1 TMP_SW_D5_20250926 TRUE
## 2 TMP_SW_H3_20250929 TRUE
## 3 TMP_SW_H6_20250926 TRUE
```

### Export Processed Data

```
## Export Processed Data
```

```
## # A tibble: 3 x 21
##   Project      plot grid Depth_cm sample_type Vial_ID date npoc_mgL npoc_uM
##   <chr>        <chr> <chr>    <dbl> <chr>      <chr>  <chr>    <dbl>    <dbl>
## 1 COMPASS: TEMP~ SW    D5        15 DOC      SW_D5_~ 2025~    15.7    1311.
## 2 COMPASS: TEMP~ SW    H3        15 DOC      SW_H3_~ 2025~    17.7    1475
## 3 COMPASS: TEMP~ SW    H6        15 DOC      SW_H6_~ 2025~    20.8    1731.
## # i 12 more variables: npoc_flag <chr>, tdn_mgL <dbl>, tdn_uM <dbl>,
## #   tdn_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>,
## #   Evacuation_date_YYYYMMDD <dbl>, Collection_Date_YYYYMMDD <dbl>,
## #   Collection_Start_Time_24hrs <dbl>, Collection_End_Time_24hrs <dbl>,
## #   EST_EDT <chr>, Volume_mL <dbl>
```

```
#end
```

## Assess Duplicates - NO DUPLICATES ON THIS RUN

```
“‘{#r Check Duplicates, echo=FALSE}
```

```
cat(“Assess Duplicates”)
```

```
#Take a look at the raw data #head(dat_raw)
```

```
#pull out any rows that have “dup” in the sample_name column dups <- dat_raw %>%
```

```
select(!c(npoc_flag, tdn_flag)) %>% filter(str_detect(sample_name, “dup”)) #have to change this to match data
```

```
#create a new dataframe and remove dups from sample dataframe dat_raw2 <- dat_raw %>%
filter(!str_detect(sample_name, “dup”))
```

```
#remove the dup from these IDs so we will have duplicate sample names dupssample_name <-
gsub(“_dup”, “”, as.character(dupssample_name)) dups <- dups[, -c(4)] #remove the run date time
for colnames(dups) <- c(“sample_name”, “npoc_raw_dup”, “tdn_raw_dup”) head(dups)
```

```
QAdups <- merge(dat_raw2, dups) head(QAdups)
```

```
df2 <- as.data.frame(QAdupsnpoc_raw)df2dups <- QAdups$npoc_raw_dup
```

```
df2sds <- apply(df2, 1, sd)df2mean <- apply(df2, 1, mean)
```

```

QAdupsnpocdupscv <- -(df2sds/df2mean)*100
QAdupsnpoc_dups_cv_flag <- ifelse(QAdups$npoc_dups_cv
<10, 'YES', 'NO, rerun')

df3 <- as.data.frame(QAdupstdnraw)
df3dups <- QAdups$tdn_raw_dup

df3sds <- apply(df3, 1, sd)
df3mean <- apply(df3, 1, mean)

QAdupstdndupscv <- -(df3sds/df3mean) * 100
QAdupstdn_dups_cv_flag <- ifelse(QAdups$tdn_dups_cv
<10, 'YES', 'NO, rerun')

head(QAdups)

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
C_dups <- ggplot(data = QAdups, aes(x = sample_name, y = npoc_dups_cv, fill = npoc_dups_cv_flag))
+ geom_bar(stat = 'identity') + theme_classic() + labs(x = "Sample ID", y = "CV of NPOC
Dups (%)") + scale_fill_manual(values = c("YES" = "darkgreen", "NO, rerun" = "red")) +
theme(legend.position = "none") + geom_hline(yintercept = 10, linetype = "dashed", color = "black", size = 1)
+ guides(fill = guide_legend(title = "CV Between Dups <10%")) + theme(axis.text.x = element_text(angle
= 90, hjust = 0.5))

N_dups <- ggplot(data = QAdups, aes(x = sample_name, y = tdn_dups_cv, fill = tdn_dups_cv_flag)) +
geom_bar(stat = 'identity') + theme_classic() + labs(x = "Sample ID", y = "CV of TN Dups (%)") +
scale_fill_manual(values = c("YES" = "darkgreen", "NO, rerun" = "red")) + theme(legend.position = "none")
+ geom_hline(yintercept = 10, linetype = "dashed", color = "black", size = 1) + guides(fill = guide_legend(title = "CV
Between Dups <10%")) + theme(axis.text.x = element_text(angle = 90, hjust = 0.5))

ggarrange(C_dups, N_dups, ncol = 2, nrow = 1)

#calculate the percent of check standards that are within the range based on the flag c_dups_percent <-
(sum(QAdupsnpocdupscv_flag == "YES")/nrow(QAdups))*100
ndupsppercent <- -(sum(QAdupstdn_dups_cv_flag
== "YES")/nrow(QAdups))*100

#report out if flags indicate need for rerun ifelse(c_dups_percent >= chks_flag, ">60% of Car-
bon Duplicates have a CV <10%", "<60% of Carbon Duplicates have a CB <10% - REASSESS")
ifelse(n_dups_percent >= chks_flag, ">60% of Nitrogen Duplicates have a CV <10%", "<60% of Nitrogen
Duplicates have a CB <10% - REASSESS")

#write out a flag to the sample dataframe if more than 60% of the dups have CVs out of range if
(c_dups_percent <= chks_flag) { dat_rawnpoc_flag <- ifelse(dat_rawnpoc_flag != "", paste0(dat_raw$npoc_flag, ";
NPOC dups out of range"), "NPOC dups out of range" ) }

if (n_dups_percent <= chks_flag) { # assuming you have tn_chks_percent similarly dat_rawtdn_flag <-
ifelse(dat_rawtdn_flag != "", paste0(dat_raw$tdn_flag, "; TN dups out of range"), "TN dups out of range"
) }

““

```