

# COMPASS: TEMPEST Discrete DOC Data QAQC

September 2024

2025-06-23

## Run Information

```
#identify which section you are in  
cat("Run Information")
```

```
## Run Information
```

```
#a link to the Gitbook or whatever protocol you are using for this analysis  
#steph will add this soon
```

```
#anything that needs to be changed do this in the first chunk
```

```
Date_Run = "09/10/24"
```

```
Run_by = "Stephanie J. Wilson"
```

```
Script_run_by = "Stephanie J. Wilson"
```

```
run_notes = "This run the samples names were not put in correctly, they are  
changed in the code to be the correct format"
```

```
#file path and name for summary file
```

```
raw_file_name = "tmp_doc_raw_data_2024/TMP_202409.txt"
```

```
#file path and name for the all peaks file
```

```
raw_allpeaks_name = "tmp_doc_raw_data_2024/TMP_202409_allpeaks.txt"
```

```
#file path and name for processed data after QAQC
```

```
processed_file_name = "tmp_doc_processed_data_2024/TMP_PW_DOC_Processed_202409.csv"
```

```
#check standard concentrations - Update if running different checks:
```

```
chk_std_c = 1
```

```
chk_std_n = 1
```

```
#Log path
```

```
Log_path = "tmp_doc_raw_data_2024/COMPASS_TMP_TOCTN_QAQClog_2024.csv"
```

## Setup

### Pull in active porewater tracking inventory sheet

```
## File already exists. No download needed.
```

## Import Data Functions

### Import Sample Data

```
## Import Sample Data
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## # A tibble: 3 x 4
```

```
##   sample_name    npoc_raw tdn_raw run_datetime
```

```
##   <chr>          <dbl>    <dbl> <chr>
```

```
## 1 TMP_SW_C3_15cm    9.33    0.666 9/10/2024 3:03:33 PM
```

```
## 2 TMP_SW_D5_15cm    5.78    0.268 9/10/2024 3:32:19 PM
```

```
## 3 TMP_SW_F6_15cm    6.01    0.353 9/10/2024 4:01:16 PM
```

## Assessing standard Curves


```
## Assess the Standard Curve
```

```
## New names:
```

```
## 'geom_smooth()' using formula = 'y ~ x'
```


```
## * '' -> '...18'
```

NPOC Std Curve by Date




```
## 'geom_smooth()' using formula = 'y ~ x'
```

TN Std Curve by Date



```
## Warning: Removed 15 rows containing missing values or values outside the scale range
## ('geom_point()').
```

```
## Warning: Removed 15 rows containing missing values or values outside the scale range
## ('geom_line()').
```



```
## [1] "NPOC Curve r2 GOOD"
```

```
## [1] "TN Curve r2 GOOD"
```


## Assess Check Standards

```
## Assess the Check Standards
```

```
## New names:
## * ' ' -> '...14'
```

```
## [1] "Carbon Check Standard RSD within Range"
```

```
## [1] "Nitrogen Check Standard RSD within Range"
```



```
## [1] "<60% of Carbon Check Standards are within range of the expected concentration - REASSESS"
```

```
## [1] ">60% of Nitrogen Check Standards are within range of the expected concentration"
```

## Assess Blanks


```
## Assess Blanks
```

```
## New names:
```

```
## * ' ' -> '...14'
```

```
## [1] ">60% of Carbon Blank concentrations are below the lower 25% quartile of samples"
```

```
## [1] ">60% of Nitrogen Blank concentrations are below the lower 25% quartile of samples"
```



```
## carbon blanks:
```

```
## [1] 0.3744933
```

```
## nitrogen blanks:
```


```
## [1] 0.08074067
```


**Assess Duplicates - NO Duplicates in this run**



## Sample Flagging

## Sample Flagging





## Visualize Data by Plot

## Visualize Data


##	Site_Code	Plot	Grid_Square	Date
## 1	TMP	SW	C3	20240906
## 2	TMP	SW	D5	20240906
## 3	TMP	SW	F6	20240906


##	Site_Code	Plot	Grid_Square	Date	sample_name	npoc_raw	tdn_raw
## 1	TMP	SW	C3	20240906	TMP_SW_C3_20240906	9.329	0.6664
## 2	TMP	SW	D5	20240906	TMP_SW_D5_20240906	5.781	0.2683
## 3	TMP	SW	F6	20240906	TMP_SW_F6_20240906	6.013	0.3530

##	run_datetime	npoc_flag	tdn_flag
## 1	9/10/2024 3:03:33 PM	NPOC checks out of range	
## 2	9/10/2024 3:32:19 PM	NPOC checks out of range blank is < 25% of sample value	
## 3	9/10/2024 4:01:16 PM	NPOC checks out of range	

## Carbon by Plot





## Convert data from mg/L to uMoles/L

### Add in/check metadata

```
## Check Sample IDs with Metadata
```

```
## # A tibble: 3 x 2
##   sample_name      metadata_recorded
##   <chr>            <lgl>
## 1 TMP_SW_C3_20240906 TRUE
## 2 TMP_SW_D5_20240906 TRUE
## 3 TMP_SW_F6_20240906 TRUE
```

### Export Processed Data

```
## Export Processed Data
```

```
## # A tibble: 3 x 21
##   Project      plot grid Depth_cm sample_type Vial_ID date npoc_mgL npoc_uM
##   <chr>        <chr> <chr>    <dbl> <chr>      <chr>  <chr>    <dbl>    <dbl>
## 1 COMPASS: TEMP~ SW    C3        15 DOC      SW_C3_~ 2024~      9.33     777.
## 2 COMPASS: TEMP~ SW    D5        15 DOC      SW_D5_~ 2024~      5.78     482.
## 3 COMPASS: TEMP~ SW    F6        15 DOC      SW_F6_~ 2024~      6.01     501.
## # i 12 more variables: npoc_flag <chr>, tdn_mgL <dbl>, tdn_uM <dbl>,
## #   tdn_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>,
## #   Evacuation_date_YYYYMMDD <dbl>, Collection_Date_YYYYMMDD <dbl>,
## #   Collection_Start_Time_24hrs <dbl>, Collection_End_Time_24hrs <dbl>,
## #   EST_EDT <chr>, Volume_mL <dbl>
```

```
#end
```