

Dionex_COMPASS_TEMPEST

Stephanie J. Wilson

2023-06-07

Daily Set up

Read in and Format the raw data - change wd & file names (Make sure files are UTF-8 encoded)

```
# SULFATE DATA:
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Sdat <- read.table("Raw Data/COMPASS_TEMPEST_20230607_S04.txt", sep='\t' , header=T, skip=3)  
head(Sdat)
```

```
##      X          X.1          X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4  
## 1 1 Lab Blank      Unknown 4.347  0.0005  1.32  0.0006  0.00  
## 2 2 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.  
## 3 3 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.  
## 4 4 Standard 1 Calibration Standard 4.357  0.4391  6.60  0.5115  2.29  
## 5 5 Standard 2 Calibration Standard 4.353  0.9724  6.67  1.1327  5.14  
## 6 6 Standard 3 Calibration Standard 4.350  1.9875  7.45  2.3152  10.10  
##      IC.S04.5  
## 1      Rd  
## 2      n.a.  
## 3      n.a.  
## 4      M  
## 5      M  
## 6      M
```

```
## Only keep the columns that we need
```

```
Sdat <- Sdat[ ,c(2,5,7)] # dont need this here  
head(Sdat)
```

```
##      X.1 IC.S04.1 IC.S04.3  
## 1 Lab Blank  0.0005  0.0006  
## 2 Lab Blank  n.a.    n.a.  
## 3 Lab Blank  n.a.    n.a.  
## 4 Standard 1  0.4391  0.5115  
## 5 Standard 2  0.9724  1.1327  
## 6 Standard 3  1.9875  2.3152
```

```
## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "S04_ppm", "S04_Area")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat$S04_Area <- as.numeric(Sdat$S04_Area)
```

```
## Warning: NAs introduced by coercion
```

```
Sdat <- as.data.frame(Sdat)
head(Sdat)
```

```
##      Sample_ID S04_ppm S04_Area
## 1 Lab Blank    0.0005    0.0006
## 2 Lab Blank         NA         NA
## 3 Lab Blank         NA         NA
## 4 Standard 1    0.4391    0.5115
## 5 Standard 2    0.9724    1.1327
## 6 Standard 3    1.9875    2.3152
```

```
#Chloride data
```

```
## Read in raw data file from Dionex - copied and saved as a txt
```

```
Cldat <- read.table("Raw Data/COMPASS_TEMPEST_20230607_Cl.txt", sep='\t' , header=T, skip=3)
head(Cldat)
```

```
##      X      X.1      X.2 IC.Cl IC.Cl.1 IC.Cl.2 IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank      Unknown 3.420  0.0186  54.19  0.0261  0.07
## 2 2 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 3 3 Lab Blank      Unknown  n.a.    n.a.    n.a.    n.a.    n.a.
## 4 4 Standard 1 Calibration Standard 3.430  5.1187  92.59  7.1765  39.52
## 5 5 Standard 2 Calibration Standard 3.427  11.2205  92.58  15.7314  85.63
## 6 6 Standard 3 Calibration Standard 3.427  20.4308  92.19  28.6445  154.92
##      IC.Cl.5
## 1      MB
## 2      n.a.
## 3      n.a.
## 4      M
## 5      M
## 6      M
```

```
## Only keep the columns that we need
```

```
Cldat <- Cldat[ ,c(2,5,7)]
head(Cldat)
```

```
##      X.1 IC.Cl.1 IC.Cl.3
## 1 Lab Blank 0.0186 0.0261
## 2 Lab Blank  n.a.    n.a.
## 3 Lab Blank  n.a.    n.a.
## 4 Standard 1 5.1187 7.1765
## 5 Standard 2 11.2205 15.7314
## 6 Standard 3 20.4308 28.6445
```

```
## Name the columns correctly
colnames(Cldat) <- c("Sample_ID", "Cl_ppm", "Cl_Area")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat$Cl_Area <- as.numeric(Cldat$Cl_Area)
```

```
## Warning: NAs introduced by coercion
```

```
Cldat <- as.data.frame(Cldat)
head(Cldat)
```

```
##      Sample_ID Cl_ppm Cl_Area
## 1 Lab Blank  0.0186  0.0261
## 2 Lab Blank      NA      NA
## 3 Lab Blank      NA      NA
## 4 Standard 1  5.1187  7.1765
## 5 Standard 2 11.2205 15.7314
## 6 Standard 3 20.4308 28.6445
```

```
## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1          557.6408  23.3103 4194.5484 280.6096
## 2          557.6408  23.3103  979.1875  68.8194
## 3          557.6408  23.3103   0.0014   0.0020
## 4          557.6408  23.3103 1429.0768  72.1471
## 5          557.6408  23.3103      NA      NA
## 6          140.4325   7.6988 4194.5484 280.6096
```

```
## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID==""), ]
head(all_dat)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 26 Lab Blank  5e-04  6e-04 0.0186  0.0261
## 27 Lab Blank  5e-04  6e-04      NA      NA
## 28 Lab Blank  5e-04  6e-04      NA      NA
## 29 Lab Blank  5e-04  6e-04 0.0014  0.0020
## 30 Lab Blank  5e-04  6e-04 0.0162  0.0227
## 31 Lab Blank  5e-04  6e-04 0.0093  0.0131
```

```
all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 147 Standard 1  0.4391  0.5115 5.1187 7.1765
## 148 Standard 1  0.4391  0.5115 5.1585 7.2323
## 149 Standard 1  0.4391  0.5115 5.1432 7.2109
## 150 Standard 1  0.4479  0.5218 5.1187 7.1765
## 151 Standard 1  0.4479  0.5218 5.1585 7.2323
## 152 Standard 1  0.4479  0.5218 5.1432 7.2109
```

```
all_dat[is.na(all_dat)] <- 0
```

Pull out standards - could do some checks here if we want

```
stds <- all_dat[grepl("Standard", all_dat$Sample_ID),]
head(stds)
```

```
##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 147 Standard 1  0.4391  0.5115 5.1187 7.1765
## 148 Standard 1  0.4391  0.5115 5.1585 7.2323
## 149 Standard 1  0.4391  0.5115 5.1432 7.2109
## 150 Standard 1  0.4479  0.5218 5.1187 7.1765
## 151 Standard 1  0.4479  0.5218 5.1585 7.2323
## 152 Standard 1  0.4479  0.5218 5.1432 7.2109
```

```
stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(S04_ppm), sd = sd(S04_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_S)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean      sd      cv flag
##   <fct>      <dbl>   <dbl>   <dbl> <chr>
## 1 Standard 1  0.452 0.0136  3.02 NO, rerun
## 2 Standard 2  0.984 0.0162  1.64 YES
## 3 Standard 3  1.94 0.0543  2.80 NO, rerun
## 4 Standard 4  8.41 1.83    21.8 NO, rerun
## 5 Standard 5 20.0 0.00254 0.0127 YES
```

```
stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv < 2, 'YES', 'NO, rerun')
head(stds_chk_Cl)
```

```
## # A tibble: 5 x 5
##   Sample_ID mean      sd      cv flag
##   <fct>      <dbl>   <dbl>   <dbl> <chr>
## 1 Standard 1  5.14 0.0174  0.338 YES
```

```
## 2 Standard 2 11.2 0.00524 0.0467 YES
## 3 Standard 3 20.4 0.00439 0.0215 YES
## 4 Standard 4 101. 0.0511 0.0506 YES
## 5 Standard 5 200. 0.407 0.204 YES
```

```
lmS <- lm(stds$S04_Area ~ stds$S04_ppm)
S04_sum <- summary(lmS)
S04_Slope <- S04_sum$coefficients[2, 1]
S04_Int <- S04_sum$coefficients[1, 1]
```

```
lmCl <- lm(stds$Cl_Area ~ stds$Cl_ppm)
Cl_sum <- summary(lmCl)
Cl_Slope <- Cl_sum$coefficients[2, 1]
Cl_Int <- Cl_sum$coefficients[1, 1]
```

Calculate mmol/L concentrations

```
#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)
```

```
##                               Sample_ID S04_ppm S04_Area Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182      TEMPEST_FW_B4_20230606      2.2280   2.0763 173.8469
##      Cl_Area
## 177 117.2776
## 178 117.6170
## 179 101.9715
## 180 110.1445
## 181 109.4211
## 182 194.9897
```

```
# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06      #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)
```

```
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000

head(sampledat)
```

```
##                               Sample_ID S04_ppm S04_Area  Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527  12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408  12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045  11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515  12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687  12.0490 3902.2588
## 182      TEMPEST_FW_B4_20230606      2.2280    2.0763  173.8469
##      Cl_Area      S04_mM      Cl_mM  salinity
## 177 117.2776 17.3441266 117.981489 7.5577019
## 178 117.6170 17.3936619 118.322945 7.5795750
## 179 101.9715 14.9939021 102.583532 6.5713373
## 180 110.1445 16.1900031 110.805585 7.0980268
## 181 109.4211 16.1312757 110.077822 7.0514077
## 182 194.9897  0.0694947   4.904003 0.3141673
```

Pull out dups and check with percent difference

```
#Show me the data that we have from the calculations
head(sampledat)
```

```
##                               Sample_ID S04_ppm S04_Area  Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527  12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408  12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045  11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515  12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687  12.0490 3902.2588
## 182      TEMPEST_FW_B4_20230606      2.2280    2.0763  173.8469
##      Cl_Area      S04_mM      Cl_mM  salinity
## 177 117.2776 17.3441266 117.981489 7.5577019
## 178 117.6170 17.3936619 118.322945 7.5795750
## 179 101.9715 14.9939021 102.583532 6.5713373
## 180 110.1445 16.1900031 110.805585 7.0980268
## 181 109.4211 16.1312757 110.077822 7.0514077
## 182 194.9897  0.0694947   4.904003 0.3141673
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "_dup")) #have to change this to match data
head(dups)
```

```
##                               Sample_ID S04_ppm S04_Area  Cl_ppm  Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687  12.0490 3902.259 109.4211
## 2      TEMPEST_SW_I5_20230606_dup 240.9080    5.6127 2335.150  65.4786
##      S04_mM      Cl_mM salinity
```

```
## 1 16.131276 110.07782 7.051408
## 2 7.514286 65.87165 4.219642
```

```
#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "_dup")) %>%
  filter(!str_detect(Sample_ID, "_spk"))
head(sampledat2)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1 TEMPEST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438 117.2776
## 2 TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484 117.6170
## 3 TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862 101.9715
## 4 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580 110.1445
## 5 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469 194.9897
## 6 TEMPEST_FW_C3_20230606 12.2001 12.9198 27.3401 34.8467
## S04_mM Cl_mM salinity
## 1 17.3441266 117.9814894 7.55770195
## 2 17.3936619 118.3229450 7.57957496
## 3 14.9939021 102.5835317 6.57133726
## 4 16.1900031 110.8055853 7.09802681
## 5 0.0694947 4.9040028 0.31416735
## 6 0.3805396 0.7712299 0.04942956
```

```
#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[, -c(2:5)]
colnames(dups) <- c('Sample_ID', 'S04_mM_dup', 'Cl_mM_dup', 'salinity_dup')
head(dups)
```

```
##           Sample_ID S04_mM_dup Cl_mM_dup salinity_dup
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 16.131276 110.07782 7.051408
## 2 TEMPEST_SW_I5_20230606 7.514286 65.87165 4.219642
```

```
#put it back together with the old data set and look for duplicates
QAdups <- merge(sampledat2, dups)
head(QAdups)
```

```
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.058 110.1445
## 2 TEMPEST_SW_I5_20230606 236.2910 5.5051 2285.850 64.0963
## S04_mM Cl_mM salinity S04_mM_dup Cl_mM_dup salinity_dup
## 1 16.190003 110.80559 7.098027 16.131276 110.07782 7.051408
## 2 7.370274 64.48097 4.130558 7.514286 65.87165 4.219642
```

```
QAdups$S04_dups_chk <- ((abs(QAdups$S04_mM-QAdups$S04_mM_dup))/((QAdups$S04_mM+QAdups$S04_mM_dup)/2))*100
QAdups$S04_dups_flag <- ifelse(QAdups$S04_dups_chk <10, 'YES', 'NO, rerun')
```

```
QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')
```

```
head(QAdups)
```

```
##                               Sample_ID  S04_ppm S04_Area  Cl_ppm  Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515  12.0928 3928.058 110.1445
## 2                               TEMPEST_SW_I5_20230606 236.2910   5.5051 2285.850  64.0963
##      S04_mM      Cl_mM salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk
## 1 16.190003 110.80559 7.098027 16.131276 110.07782    7.051408    0.3633977
## 2  7.370274  64.48097 4.130558  7.514286  65.87165    4.219642    1.9350418
##      S04_dups_flag Cl_dups_chk Cl_dups_flag
## 1                YES    0.6589567                YES
## 2                YES    2.1337241                YES
```


#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this

```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)
```

```
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```

```
Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)
```

```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES         2  YES         2      2        100        100
```

Pull out dups and check with cv

```
#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100
```

```
df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2,1,sd)

QAdups$S04_dups_cv <- (df2$sds)/((QAdups$S04_mM+QAdups$S04_mM_dup)/2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv <11, 'YES', 'NO, rerun')

head(QAdups)
```


```
##                               Sample_ID  S04_ppm S04_Area  Cl_ppm  Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515  12.0928 3928.058 110.1445
## 2                               TEMPEST_SW_I5_20230606 236.2910   5.5051 2285.850  64.0963
##      S04_mM      Cl_mM salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk
## 1 16.190003 110.80559 7.098027  16.131276 110.07782    7.051408    0.3633977
## 2  7.370274  64.48097 4.130558   7.514286  65.87165    4.219642    1.9350418
##      S04_dups_flag Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1              YES  0.6589567            YES  0.256961            YES
## 2              YES  2.1337241            YES  1.368281            YES
```

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this

```
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
  color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar,ncol=2, nrow=1)
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES          2  YES          2      2         100       100
```

Pull out spikes and check with dionex output conc.

```
##Merge sample data with dilution data
Dilutions <- read.csv("Raw Data/TMP 2023 Sample Dilutions from Dionex.csv")
```

```
sampldat_dils <- merge(sampldat, Dilutions, by = "Sample_ID")
unique(sampldat_dils$Dilution)
```

```
## [1] 50.00 1.25 1.10
```

```
#Show me the data that we have from the calculations
head(sampldat_dils)
```

```
##
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1 TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438 117.2776
## 2 TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484 117.6170
## 3 TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862 101.9715
## 4 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580 110.1445
## 5 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588 109.4211
## 6 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469 194.9897
##      S04_mM      Cl_mM salinity Dilution
## 1 17.3441266 117.981489 7.5577019 50.00
## 2 17.3936619 118.322945 7.5795750 50.00
## 3 14.9939021 102.583532 6.5713373 50.00
## 4 16.1900031 110.805585 7.0980268 50.00
## 5 16.1312757 110.077822 7.0514077 50.00
## 6 0.0694947 4.904003 0.3141673 1.25
```

```
#pull out any rows that have "spk" in the SampleID column
spks <- sampldat_dils %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)
```

```
##
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1 TEMPEST_FW_I5_20230606_spk 16.5901 15.4606 31.4901 35.3198
## 2 TEMPEST_FW_SOURCE_20230606_1511_spk 16.3183 17.2810 24.0314 30.6296
##      S04_mM      Cl_mM salinity Dilution
## 1 0.5174704 0.8882962 0.05692861 1.25
## 2 0.5089925 0.6778956 0.04345074 1.10
```

```
#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[,c(1,6)]
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
head(spks)
```

```
##
##           Sample_ID S04_mM_spk
## 1 TEMPEST_FW_I5_20230606 0.5174704
## 2 TEMPEST_FW_SOURCE_20230606_1511 0.5089925
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampldat_dils, spks)
head(QAspks)
```

```
##
##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area S04_mM
```

```
## 1          TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##          Cl_mM    salinity Dilution SO4_mM_spk
## 1 0.8921834 0.05717762      1.25 0.5174704
## 2 0.6922793 0.04437213      1.10 0.5089925
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10              # in uL
spkvol <- spkvol/1000000
#spike for these samples was 10uL of the 250mM standard
QAspks$SO4_spk_Conc <- (spkconc)*spkvol      # mmoles of SO4
head(QAspks)
```

```
##          Sample_ID SO4_ppm SO4_Area Cl_ppm Cl_Area    SO4_mM
## 1          TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##          Cl_mM    salinity Dilution SO4_mM_spk SO4_spk_Conc
## 1 0.8921834 0.05717762      1.25 0.5174704 7.797879e-05
## 2 0.6922793 0.04437213      1.10 0.5089925 7.797879e-05
```

```
#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
# QAspks$Dilution <- 1
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
# QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
unique(QAspks$Dilution)
```

```
## [1] 1.25 1.10
```

```
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(QAspks$Dilution == 50, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.85, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 51.16, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 49.27, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.51, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 5.15, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 1.25, 1500, QAspks$SampleVol)
```

```

QAspks$SampleVol <- ifelse(QAspks$Dilution == 1.1, 1500, QAspks$SampleVol)

# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 1483, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1      TEMPEST_FW_I5_20230606 12.8534  11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474  14.7702 24.5413 31.2795 0.4350405
##           Cl_mM  salinity Dilution S04_mM_spk S04_spk_Conc SampleVol
## 1 0.8921834 0.05717762      1.25  0.5174704 7.797879e-05      0.0015
## 2 0.6922793 0.04437213      1.10  0.5089925 7.797879e-05      0.0015

#gives us the total S04 in the sample in mmoles
QAspks$S04_Total_unspkd <- (QAspks$S04_mM/QAspks$Dilution)*(QAspks$SampleVol)

##total S04 in spiked sample in mmoles
QAspks$S04_Total_spkd <- (QAspks$S04_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)

QAspks$S04_expctd_spkd <- (QAspks$S04_Total_unspkd + QAspks$S04_spk_Conc)
QAspks$spk_recovery <- (QAspks$S04_Total_spkd/QAspks$S04_expctd_spkd)*100
QAspks$S04_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun')

head(QAspks)


##           Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1      TEMPEST_FW_I5_20230606 12.8534  11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474  14.7702 24.5413 31.2795 0.4350405
##           Cl_mM  salinity Dilution S04_mM_spk S04_spk_Conc SampleVol
## 1 0.8921834 0.05717762      1.25  0.5174704 7.797879e-05      0.0015
## 2 0.6922793 0.04437213      1.10  0.5089925 7.797879e-05      0.0015
##   S04_Total_unspkd S04_Total_spkd S04_expctd_spkd spk_recovery S04_spks_flag
## 1      0.0004811004  0.0006251042    0.0005590792     111.8096         YES
## 2      0.0005932371  0.0006987079    0.0006712159     104.0959         YES

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +

```

```
scale_fill_manual(values=c("lightblue4")) +
theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
theme(legend.position="none") +
geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

spksbar



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 YES           2      2    100
```

Pull out spikes and check with area calc

```
##Merge sample data with dilution data
```

```
Dilutions <- read.csv("Raw Data/TMP 2023 Sample Dilutions from Dionex.csv")
sampledat_dils <- merge(sampledat, Dilutions, by = "Sample_ID")
unique(sampledat_dils$Dilution)
```

```
## [1] 50.00 1.25 1.10
```

```
#Show me the data that we have from the calculations
```

```
head(sampledat_dils)
```

```
##           Sample_ID S04_ppm S04_Area   Cl_ppm Cl_Area
## 1 TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438 117.2776
## 2 TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484 117.6170
## 3 TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862 101.9715
## 4 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580 110.1445
## 5 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588 109.4211
## 6 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469 194.9897
##      S04_mM      Cl_mM salinity Dilution
## 1 17.3441266 117.981489 7.5577019 50.00
## 2 17.3936619 118.322945 7.5795750 50.00
## 3 14.9939021 102.583532 6.5713373 50.00
## 4 16.1900031 110.805585 7.0980268 50.00
## 5 16.1312757 110.077822 7.0514077 50.00
## 6 0.0694947 4.904003 0.3141673 1.25
```

```
#now use area of all samples and calculate undilution corrected concentrations in ug/mL
```

```
sampledat_dils$S04_ugmL <- ((sampledat_dils$S04_Area)-S04_Int)/S04_Slope
sampledat_dils$Cl_ugmL <- (sampledat_dils$Cl_Area-Cl_Int)/Cl_Slope
head(sampledat_dils)
```

```
##           Sample_ID S04_ppm S04_Area   Cl_ppm Cl_Area
## 1 TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438 117.2776
## 2 TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484 117.6170
## 3 TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862 101.9715
## 4 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580 110.1445
## 5 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588 109.4211
## 6 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469 194.9897
##      S04_mM      Cl_mM salinity Dilution S04_ugmL Cl_ugmL
## 1 17.3441266 117.981489 7.5577019 50.00 11.121064 83.64890
## 2 17.3936619 118.322945 7.5795750 50.00 11.152826 83.89098
## 3 14.9939021 102.583532 6.5713373 50.00 9.614066 72.73174
## 4 16.1900031 110.805585 7.0980268 50.00 10.381000 78.56118
## 5 16.1312757 110.077822 7.0514077 50.00 10.343400 78.04521
## 6 0.0694947 4.904003 0.3141673 1.25 1.782401 139.07750
```

```
#pull out any rows that have "spk" in the SampleID column
```

```
spks <- sampledat_dils %>%
  filter(str_detect(Sample_ID, "spk")) #have to change this to match data
head(spks)
```



```
##                               Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1          TEMPEST_FW_I5_20230606_spk 16.5901  15.4606 31.4901 35.3198
## 2 TEMPEST_FW_SOURCE_20230606_1511_spk 16.3183  17.2810 24.0314 30.6296
##      S04_mM      Cl_mM      salinity Dilution S04_u/mL  Cl_u/mL
## 1 0.5174704 0.8882962 0.05692861      1.25 13.27207 25.19204
## 2 0.5089925 0.6778956 0.04345074      1.10 14.83478 21.84672
```

```
## Only keep the columns that we need
spks <- spks[,c(1,9,10)]
head(spks)
```

```
##                               Sample_ID Dilution S04_u/mL
## 1          TEMPEST_FW_I5_20230606_spk      1.25 13.27207
## 2 TEMPEST_FW_SOURCE_20230606_1511_spk      1.10 14.83478
```

```
#remove the spk from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
colnames(spks) <- c('Sample_ID', 'Dilution', 'S04_u/mL_spk')
head(spks)
```

```
##                               Sample_ID Dilution S04_u/mL_spk
## 1          TEMPEST_FW_I5_20230606      1.25      13.27207
## 2 TEMPEST_FW_SOURCE_20230606_1511      1.10      14.83478
```

```
#put it back together with the old data set and look for duplicates
QAspks <- merge(sampldat_dils, spks)
head(QAspks)
```

```
##                               Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area
## 1          TEMPEST_FW_I5_20230606      1.25 12.8534  11.9783 31.6279 35.4744
## 2 TEMPEST_FW_SOURCE_20230606_1511      1.10 13.9474  14.7702 24.5413 31.2795
##      S04_mM      Cl_mM      salinity S04_u/mL  Cl_u/mL S04_u/mL_spk
## 1 0.4009170 0.8921834 0.05717762 10.28271 25.30230      13.27207
## 2 0.4350405 0.6922793 0.04437213 12.67940 22.31027      14.83478
```

```
#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250)      # in ug
spkvol <- 10          # in uL
spkvol <- spkvol/1000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol      # mmoles of S04
head(QAspks)
```

```
##                               Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area
## 1          TEMPEST_FW_I5_20230606      1.25 12.8534  11.9783 31.6279 35.4744
## 2 TEMPEST_FW_SOURCE_20230606_1511      1.10 13.9474  14.7702 24.5413 31.2795
##      S04_mM      Cl_mM      salinity S04_u/mL  Cl_u/mL S04_u/mL_spk S04_spk_Conc
## 1 0.4009170 0.8921834 0.05717762 10.28271 25.30230      13.27207      2.5
## 2 0.4350405 0.6922793 0.04437213 12.67940 22.31027      14.83478      2.5
```

```

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(QAspks$Dilution == 50, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 1.1, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 5.15, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 4.53, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 56.51, 1483, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 1.25, 1500, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(QAspks$Dilution == 1.1, 1500, QAspks$SampleVol)

# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "TEMPEST"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
# QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000
head(QAspks)

```

```

##           Sample_ID Dilution S04_ppm S04_Area Cl_ppm Cl_Area
## 1      TEMPEST_FW_I5_20230606      1.25 12.8534  11.9783 31.6279 35.4744
## 2 TEMPEST_FW_SOURCE_20230606_1511      1.10 13.9474  14.7702 24.5413 31.2795
##      S04_mM      Cl_mM      salinity S04_ugmL Cl_ugmL S04_ugmL_spk S04_spk_Conc
## 1 0.4009170 0.8921834 0.05717762 10.28271 25.30230      13.27207      2.5

```

```
## 2 0.4350405 0.6922793 0.04437213 12.67940 22.31027      14.83478      2.5
##   SampleVol
## 1      1.5
## 2      1.5
```

```
#gives us the total SO4 in the sample in mmoles
```

```
QAspks$SO4_Total_unspkd <- QAspks$SO4_ugmL*QAspks$SampleVol
```

```
##total SO4 in spiked sample in mmoles
```

```
QAspks$SO4_Total_spkd <- (QAspks$SO4_ugmL_spk)*(QAspks$SampleVol+spkvol)
```

```
QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)
```

```
QAspks$spk_recovery <- (QAspks$SO4_Total_spkd/QAspks$SO4_expctd_spkd)*100
```

```
QAspks$SO4_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun
```


```
head(QAspks)
```

```
##           Sample_ID Dilution SO4_ppm SO4_Area Cl_ppm Cl_Area
## 1 TEMPEST_FW_I5_20230606      1.25 12.8534  11.9783 31.6279 35.4744
## 2 TEMPEST_FW_SOURCE_20230606_1511      1.10 13.9474  14.7702 24.5413 31.2795
##      SO4_mM      Cl_mM      salinity SO4_ugmL Cl_ugmL SO4_ugmL_spk SO4_spk_Conc
## 1 0.4009170 0.8921834 0.05717762 10.28271 25.30230      13.27207      2.5
## 2 0.4350405 0.6922793 0.04437213 12.67940 22.31027      14.83478      2.5
##   SampleVol SO4_Total_unspkd SO4_Total_spkd SO4_expctd_spkd spk_recovery
## 1      1.5      15.42406      20.04082      17.92406      111.8096
## 2      1.5      19.01909      22.40051      21.51909      104.0960
##   SO4_spks_flag
## 1      YES
## 2      YES
```

```
#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work o
```

```
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)
```

```
spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>         <int> <int>   <dbl>
## 1 YES                2     2    100
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_ ', fixed=TRUE)))
```

```
## Warning in rbind(c("TEMPEST", "SW", "SOURCE", "20230606", "1301"), c("TEMPEST",
## : number of columns of result is not a multiple of vector length (arg 5)
```

```
colnames(IDs) <- c("Project" , "Plot","Grid", "Date", "Time")
head(IDs)
```

```
##   Project   Plot   Grid   Date   Time
## 1 TEMEPST    SW SOURCE 20230606 1301
## 2 TEMPEST ESTUARY SOURCE 20230606 0632
## 3 TEMPEST ESTUARY SOURCE 20230606 1236
## 4 TEMPEST ESTUARY SOURCE 20230606 1524
## 5 TEMPEST    FW    B4 20230606 TEMPEST
## 6 TEMPEST    FW    C3 20230606 TEMPEST
```

```
#rejoin them to the dataframe
alldat <- cbind(IDs, sampled2)
head(alldat)
```

```
##   Project   Plot   Grid   Date   Time   Sample_ID
## 1 TEMEPST    SW SOURCE 20230606 1301   TEMEPST_SW_SOURCE_20230606_1301
## 2 TEMPEST ESTUARY SOURCE 20230606 0632 TEMPEST_ESTUARY_SOURCE_20230606_0632
## 3 TEMPEST ESTUARY SOURCE 20230606 1236 TEMPEST_ESTUARY_SOURCE_20230606_1236
## 4 TEMPEST ESTUARY SOURCE 20230606 1524 TEMPEST_ESTUARY_SOURCE_20230606_1524
## 5 TEMPEST    FW    B4 20230606 TEMPEST   TEMPEST_FW_B4_20230606
## 6 TEMPEST    FW    C3 20230606 TEMPEST   TEMPEST_FW_C3_20230606
##   SO4_ppm SO4_Area   Cl_ppm Cl_Area   SO4_mM   Cl_mM   salinity
## 1 556.0527 12.9549 4182.4438 117.2776 17.3441266 117.9814894 7.55770195
## 2 557.6408 12.9919 4194.5484 117.6170 17.3936619 118.3229450 7.57957496
## 3 480.7045 11.1994 3636.5862 101.9715 14.9939021 102.5835317 6.57133726
## 4 519.0515 12.0928 3928.0580 110.1445 16.1900031 110.8055853 7.09802681
## 5   2.2280   2.0763 173.8469 194.9897   0.0694947   4.9040028 0.31416735
## 6 12.2001 12.9198   27.3401  34.8467   0.3805396   0.7712299 0.04942956
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_TEMPEST_Dionex_ClSo4_Final_20230607.csv")
```

#Cha

END