COMPASS: TEMPEST Discrete DOC Data QAQC

July 2025

2025-10-02

Run Information

```
#identify which section you are in
cat("Run Information")
```

Run Information

```
#a link to the Gitbook or whatever protocol you are using for this analysis
  #steph will add this soon
#anything that needs to be changed do this in the first chunk
 Date Run = "07/03/25"
  Run_by = "Stephanie J. Wilson"
  Script_run_by = "Stephanie J. Wilson"
 run_notes = " "
  #file path and name for summary file
   raw_file_name = "tmp_doc_raw_data_2025/TMP_202507.txt"
  #file path and name for the all peaks file
   raw_allpeaks_name = "tmp_doc_raw_data_2025/TMP_202507_allpeaks.txt"
  #file path and name for processed data after QAQC
   processed_file_name = "tmp_doc_processed_data_2025/TMP_PW_DOC_Processed_202507.csv"
#check standard concentrations - Update if running different checks:
   chk std c = 50
  chk_std_n = 2
#Log path
   Log_path = "tmp_doc_raw_data_2025/COMPASS_TMP_TOCTN_QAQClog_2025.csv"
```

Setup

Pull in active porewater tracking inventory sheet

File already exists. No download needed.

Import Data Functions

Import Sample Data


```
## Import Sample Data
## New names:
## * '' -> '...14'
## # A tibble: 6 x 4
   sample_name
                       npoc_raw tdn_raw run_datetime
##
    <chr>
                        <dbl> <dbl> <chr>
## 1 TMP_C_F6_20250711
                          24.8 0.824 7/14/2025 9:35:40 PM
## 2 TMP_C_H3_20250711
                          25.4 0.926 7/14/2025 10:04:08 PM
## 3 TMP_C_H3_20250711_dup
                          25.7 1.01 7/14/2025 10:33:15 PM
                           25.8 0.862 7/14/2025 11:45:04 PM
## 4 TMP_C_H6_20250711
                      18.0 0.584 7/15/2025 12:12:38 AM
## 5 TMP_C_I5_20250711
## 6 TMP_SW_B4_20250711
                          75.8 2.60 7/15/2025 12:34:38 AM
```

Assessing standard Curves

Assess the Standard Curve


```
## New names:
## 'geom_smooth()' using formula = 'y ~ x'
## * '' -> '...18'
```

NPOC Std Curve by Date




'geom_smooth()' using formula = 'y ~ x'

TN Std Curve by Date



Slope Drift Assessment



- ## [1] "NPOC Curve r2 GOOD"
- ## [1] "TN Curve r2 GOOD"


Assess Check Standards

Assess the Check Standards

New names: ## * '' -> '...14'

[1] "Carbon Check Standard RSD within Range"

[1] "Nitrogen Check Standard RSD within Range"



[1] ">60% of Carbon Check Standards are within range of the expected concentration"

[1] ">60% of Nitrogen Check Standards are within range of the expected concentration"

Assess Blanks


Assess Blanks

New names:

* '' -> '...14'

[1] ">60% of Carbon Blank concentrations are below the lower 25% quartile of samples"

[1] ">60% of Nitrogen Blank concentrations are below the lower 25% quartile of samples"



carbon blanks:

[1] 0.5250037


nitrogen blanks:

[1] 0.03449812

Assess Duplicates - if there are any

Assess Duplicates

```
## # A tibble: 3 x 3
##
    sample_name
                       npoc_raw_dup tdn_raw_dup
    <chr>
                          <dbl>
##
                                          <dbl>
## 1 TMP C H3 20250711
                              25.7
                                          1.01
## 2 TMP_SW_D5_20250711
                                          1.27
                               31.9
## 3 TMP_FW_C6_20250711
                               10.2
                                          0.437
##
           sample_name npoc_raw tdn_raw
                                                run_datetime npoc_flag tdn_flag
## 1 TMP_C_H3_20250711
                          25.42 0.9261 7/14/2025 10:04:08 PM
                          10.30 0.4455 7/15/2025 5:55:18 AM
## 2 TMP_FW_C6_20250711
## 3 TMP_SW_D5_20250711
                          32.21 1.2600 7/15/2025 1:22:53 AM
## npoc_raw_dup tdn_raw_dup
## 1
           25.69
                      1.0120
## 2
           10.19
                      0.4368
## 3
           31.88
                      1.2750
           sample_name npoc_raw tdn_raw
                                        run_datetime npoc_flag tdn_flag
                          25.42 0.9261 7/14/2025 10:04:08 PM
## 1 TMP_C_H3_20250711
## 2 TMP_FW_C6_20250711
                          10.30 0.4455 7/15/2025 5:55:18 AM
## 3 TMP_SW_D5_20250711
                          32.21 1.2600 7/15/2025 1:22:53 AM
## npoc_raw_dup tdn_raw_dup npoc_dups_cv npoc_dups_cv_flag tdn_dups_cv
## 1
           25.69
                      1.0120
                                 1.116464
                                                       YES
                                                              9.116356
## 2
           10.19
                      0.4368
                                 1.134518
                                                       YES
                                                              2.077264
                                1.088307
## 3
           31.88
                      1.2750
                                                       YES
                                                              1.249989
## tdn_dups_cv_flag
## 1
                 YES
## 2
                 YES
## 3
                 YES
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```





[1] ">60% of Carbon Duplicates have a CV <10%"

[1] ">60% of Nitrogen Duplicates have a CV <10%"

Sample Flagging


Sample Flagging






Visualize Data by Plot

Visualize Data ## Warning in rbind(c("TMP", "C", "F6", "20250711"), c("TMP", "C", "H3", ## "20250711": number of columns of result is not a multiple of vector length (arg ## 1) Site_Code Plot Grid_Square Date NA ## 1 TMP F6 20250711 TMP C ## 2 TMP С H3 20250711 TMP ## 3 TMP C H3 20250711 dup ## 4 TMP C H6 20250711 TMP ## 5 TMP C I5 20250711 TMP B4 20250711 TMP ## 6 TMP SW Site_Code Plot Grid_Square Date NA sample_name npoc_raw TMP_C_F6_20250711 ## 1 TMP F6 20250711 TMP 24.78 C ## 2 TMP С H3 20250711 TMP TMP C H3 20250711 25.42 ## 3 TMP C H3 20250711 dup TMP_C_H3_20250711_dup 25.69 ## 4 TMP C H6 20250711 TMP TMP_C_H6_20250711 25.83 ## 5 С TMP_C_I5_20250711 TMP I5 20250711 TMP 17.98 ## 6 TMP SW B4 20250711 TMP TMP_SW_B4_20250711 75.80 run_datetime npoc_flag tdn_flag ## ${ t dn}_{ t raw}$ ## 1 0.8238 7/14/2025 9:35:40 PM ## 2 0.9261 7/14/2025 10:04:08 PM ## 3 1.0120 7/14/2025 10:33:15 PM ## 4 0.8619 7/14/2025 11:45:04 PM ## 5 0.5845 7/15/2025 12:12:38 AM ## 6 2.5970 7/15/2025 12:34:38 AM



Nitrogen by Plot



Convert data from mg/L to uMoles/L

Add in/check metadata

Check Sample IDs with Metadata

```
## # A tibble: 17 x 2
##
      sample_name
                         metadata_recorded
##
      <chr>
                         <1g1>
##
   1 TMP_C_F6_20250711
                         TRUE
##
   2 TMP C H3 20250711
                         TRUE
##
   3 TMP_C_H6_20250711
                         TRUE
  4 TMP C I5 20250711
                         TRUE
## 5 TMP_SW_B4_20250711 TRUE
## 6 TMP_SW_C3_20250711 TRUE
##
  7 TMP_SW_D5_20250711 TRUE
## 8 TMP SW E3 20250711 TRUE
## 9 TMP_SW_F6_20250711 TRUE
## 10 TMP_SW_H3_20250711 TRUE
## 11 TMP_SW_H6_20250711 TRUE
## 12 TMP_SW_I5_20250711 TRUE
## 13 TMP FW C3 20250711 TRUE
## 14 TMP FW C6 20250711 TRUE
## 15 TMP_FW_D5_20250711 TRUE
## 16 TMP_FW_F4_20250711 TRUE
## 17 TMP_FW_I5_20250711 TRUE
```

Export Processed Data

Export Processed Data

```
## # A tibble: 6 x 21
     Project
                    plot grid Depth_cm sample_type Vial_ID date npoc_mgL npoc_uM
##
     <chr>
                    <chr> <chr>
                                    <dbl> <chr>
                                                      <chr>
                                                               <chr>>
                                                                        <dbl>
                                                                                <dbl>
## 1 COMPASS: TEMP~ C
                          F6
                                       15 DOC
                                                      C_F6_D~ 2025~
                                                                         24.8
                                                                                2065
## 2 COMPASS: TEMP~ C
                                       15 DOC
                                                      C_H3_D~ 2025~
                                                                         25.4
                          НЗ
                                                                                2118.
## 3 COMPASS: TEMP~ C
                          Н6
                                       15 DOC
                                                      C H6 D~ 2025~
                                                                         25.8
                                                                                2152.
                                                      C_I5_D~ 2025~
## 4 COMPASS: TEMP~ C
                          I5
                                       15 DOC
                                                                         18.0
                                                                                1498.
## 5 COMPASS: TEMP~ SW
                          В4
                                       15 DOC
                                                      SW_B4_~ 2025~
                                                                         75.8
                                                                                6317.
## 6 COMPASS: TEMP~ SW
                          СЗ
                                       15 DOC
                                                      SW_C3_~ 2025~
                                                                         96.4
                                                                                8030
## # i 12 more variables: npoc_flag <chr>, tdn_mgL <dbl>, tdn_uM <dbl>,
       tdn_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>,
## #
       Evacuation date YYYMMDD <dbl>, Collection Date YYYYMMDD <dbl>,
## #
       Collection_Start_Time_24hrs <dbl>, Collection_End_Time_24hrs <dbl>,
## #
## #
       EST_EDT <chr>, Volume_mL <dbl>
```

#end