

Dionex_COMPASS_TEMPEST

Stephanie J. Wilson

2023-06-07

Daily Set up

Read in and Format the raw data - change wd & file names (Make sure files are UTF-8 encoded)

```
# SULFATE DATA:

## Read in raw data file from Dionex - copied and saved as a txt
Sdat <- read.table("Raw Data/COMPASS_TEMPEST_20230607_S04.txt",sep='\t' , header=T, skip=3)
head(Sdat)

##      X          X.1          X.2 IC.S04 IC.S04.1 IC.S04.2 IC.S04.3 IC.S04.4
## 1 1 Lab Blank Unknown 4.347  0.0005    1.32  0.0006    0.00
## 2 2 Lab Blank Unknown n.a.     n.a.     n.a.     n.a.     n.a.
## 3 3 Lab Blank Unknown n.a.     n.a.     n.a.     n.a.     n.a.
## 4 4 Standard 1 Calibration Standard 4.357  0.4391    6.60  0.5115    2.29
## 5 5 Standard 2 Calibration Standard 4.353  0.9724    6.67  1.1327    5.14
## 6 6 Standard 3 Calibration Standard 4.350  1.9875    7.45  2.3152   10.10
##      IC.S04.5
## 1      Rd
## 2      n.a.
## 3      n.a.
## 4      M
## 5      M
## 6      M

## Only keep the columns that we need
Sdat <- Sdat[,c(2,5,7)] # dont need this here
head(Sdat)

##          X.1 IC.S04.1 IC.S04.3
## 1  Lab Blank  0.0005  0.0006
## 2  Lab Blank    n.a.    n.a.
## 3  Lab Blank    n.a.    n.a.
## 4 Standard 1  0.4391  0.5115
## 5 Standard 2  0.9724  1.1327
## 6 Standard 3  1.9875  2.3152
```

```

## Name the columns correctly
colnames(Sdat) <- c( "Sample_ID", "S04_ppm", "S04_Area")
Sdat$Sample_ID <- as.factor(Sdat$Sample_ID)
Sdat$S04_ppm <- as.numeric(Sdat$S04_ppm)

## Warning: NAs introduced by coercion

Sdat$S04_Area <- as.numeric(Sdat$S04_Area)

## Warning: NAs introduced by coercion

Sdat <- as.data.frame(Sdat)
head(Sdat)

##      Sample_ID S04_ppm S04_Area
## 1    Lab Blank   0.0005   0.0006
## 2    Lab Blank      NA       NA
## 3    Lab Blank      NA       NA
## 4 Standard 1    0.4391   0.5115
## 5 Standard 2    0.9724   1.1327
## 6 Standard 3    1.9875   2.3152

#Chloride data
## Read in raw data file from Dionex - copied and saved as a txt
Cldat <- read.table("Raw Data/COMPASS_TEMPEST_20230607_Cl.txt",sep='\t' , header=T, skip=3)
head(Cldat)

##      X          X.1          X.2  IC.Cl  IC.Cl.1 IC.Cl.2  IC.Cl.3 IC.Cl.4
## 1 1 Lab Blank Unknown 3.420  0.0186   54.19  0.0261   0.07
## 2 2 Lab Blank Unknown n.a.     n.a.     n.a.     n.a.     n.a.
## 3 3 Lab Blank Unknown n.a.     n.a.     n.a.     n.a.     n.a.
## 4 4 Standard 1 Calibration Standard 3.430  5.1187   92.59  7.1765  39.52
## 5 5 Standard 2 Calibration Standard 3.427 11.2205   92.58 15.7314  85.63
## 6 6 Standard 3 Calibration Standard 3.427 20.4308   92.19 28.6445 154.92
##      IC.Cl.5
## 1      MB
## 2      n.a.
## 3      n.a.
## 4      M
## 5      M
## 6      M

## Only keep the columns that we need
Cldat <- Cldat[,c(2,5,7)]
head(Cldat)

##          X.1  IC.Cl.1  IC.Cl.3
## 1  Lab Blank  0.0186  0.0261
## 2  Lab Blank      n.a.     n.a.
## 3  Lab Blank      n.a.     n.a.
## 4 Standard 1  5.1187  7.1765
## 5 Standard 2 11.2205 15.7314
## 6 Standard 3 20.4308 28.6445

```

```

## Name the columns correctly
colnames(Cldat) <- c( "Sample_ID", "Cl_ppm", "Cl_Area")
Cldat$Sample_ID <- as.factor(Cldat$Sample_ID)
Cldat$Cl_ppm <- as.numeric(Cldat$Cl_ppm)

```

```

## Warning: NAs introduced by coercion

```

```

Cldat$Cl_Area <- as.numeric(Cldat$Cl_Area)

```

```

## Warning: NAs introduced by coercion

```

```

Cldat <- as.data.frame(Cldat)
head(Cldat)

```

```

##   Sample_ID  Cl_ppm Cl_Area
## 1 Lab Blank  0.0186  0.0261
## 2 Lab Blank      NA      NA
## 3 Lab Blank      NA      NA
## 4 Standard 1  5.1187  7.1765
## 5 Standard 2 11.2205 15.7314
## 6 Standard 3 20.4308 28.6445

```

```

## Bring the data back together:
all_dat <- merge(Sdat, Cldat, by="Sample_ID")
head(all_dat)

```

```

##   Sample_ID  S04_ppm S04_Area    Cl_ppm  Cl_Area
## 1           557.6408  23.3103 4194.5484 280.6096
## 2           557.6408  23.3103  979.1875  68.8194
## 3           557.6408  23.3103     0.0014  0.0020
## 4           557.6408  23.3103 1429.0768  72.1471
## 5           557.6408  23.3103       NA       NA
## 6          140.4325    7.6988 4194.5484 280.6096

```

```

## Remove empty lines
all_dat <- all_dat[!(is.na(all_dat$Sample_ID) | all_dat$Sample_ID=="") , ]
head(all_dat)

```

```

##   Sample_ID  S04_ppm S04_Area Cl_ppm  Cl_Area
## 26 Lab Blank  5e-04   6e-04  0.0186  0.0261
## 27 Lab Blank  5e-04   6e-04      NA      NA
## 28 Lab Blank  5e-04   6e-04      NA      NA
## 29 Lab Blank  5e-04   6e-04  0.0014  0.0020
## 30 Lab Blank  5e-04   6e-04  0.0162  0.0227
## 31 Lab Blank  5e-04   6e-04  0.0093  0.0131

```

```

all_dat <- droplevels(all_dat[!all_dat$Sample_ID == 'Lab Blank',])
head(all_dat)

```

```

##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 147 Standard 1  0.4391   0.5115 5.1187  7.1765
## 148 Standard 1  0.4391   0.5115 5.1585  7.2323
## 149 Standard 1  0.4391   0.5115 5.1432  7.2109
## 150 Standard 1  0.4479   0.5218 5.1187  7.1765
## 151 Standard 1  0.4479   0.5218 5.1585  7.2323
## 152 Standard 1  0.4479   0.5218 5.1432  7.2109

```

```
all_dat[is.na(all_dat)] <- 0
```

Pull out standards - could do some checks here if we want

```

stds <- all_dat[grep("Standard", all_dat$Sample_ID),]
head(stds)

```

```

##      Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 147 Standard 1  0.4391   0.5115 5.1187  7.1765
## 148 Standard 1  0.4391   0.5115 5.1585  7.2323
## 149 Standard 1  0.4391   0.5115 5.1432  7.2109
## 150 Standard 1  0.4479   0.5218 5.1187  7.1765
## 151 Standard 1  0.4479   0.5218 5.1585  7.2323
## 152 Standard 1  0.4479   0.5218 5.1432  7.2109

```

```

stds_chk_S <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(S04_ppm), sd = sd(S04_ppm))
stds_chk_S$cv <- (stds_chk_S$sd/stds_chk_S$mean)*100
stds_chk_S$flag <- ifelse(stds_chk_S$cv <2, 'YES', 'NO, rerun')
head(stds_chk_S)

```

```

## # A tibble: 5 x 5
##   Sample_ID     mean       sd       cv flag
##   <fct>        <dbl>     <dbl>     <dbl> <chr>
## 1 Standard 1  0.452    0.0136    3.02   NO, rerun
## 2 Standard 2  0.984    0.0162    1.64    YES
## 3 Standard 3  1.94     0.0543    2.80    NO, rerun
## 4 Standard 4  8.41     1.83     21.8    NO, rerun
## 5 Standard 5 20.0     0.00254   0.0127  YES

```

```

stds_chk_Cl <- stds %>%
  group_by(Sample_ID) %>%
  summarise(mean = mean(Cl_ppm), sd = sd(Cl_ppm))
stds_chk_Cl$cv <- (stds_chk_Cl$sd/stds_chk_Cl$mean)*100
stds_chk_Cl$flag <- ifelse(stds_chk_Cl$cv <2, 'YES', 'NO, rerun')
head(stds_chk_Cl)

```

```

## # A tibble: 5 x 5
##   Sample_ID     mean       sd       cv flag
##   <fct>        <dbl>     <dbl>     <dbl> <chr>
## 1 Standard 1  5.14    0.0174    0.338   YES

```

```

## 2 Standard 2 11.2 0.00524 0.0467 YES
## 3 Standard 3 20.4 0.00439 0.0215 YES
## 4 Standard 4 101. 0.0511 0.0506 YES
## 5 Standard 5 200. 0.407 0.204 YES

```

```

lmS <- lm(stds$S04_Area ~ stds$S04_ppm)
S04_sum <- summary(lmS)
S04_Slope <- S04_sum$coefficients[2, 1]
S04_Int <- S04_sum$coefficients[1, 1]

lmCl <- lm(stds$Cl_Area ~ stds$Cl_ppm)
Cl_sum <- summary(lmCl)
Cl_Slope <- Cl_sum$coefficients[2, 1]
Cl_Int <- Cl_sum$coefficients[1, 1]

```

Calculate mmol/L concentrations

```

#remove standards from sample dataframe
sampledat <- all_dat[!grepl("Standard", all_dat$Sample_ID),]
head(sampledat)

##                                     Sample_ID  S04_ppm S04_Area    Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182                  TEMPEST_FW_B4_20230606   2.2280  2.0763 173.8469
##          Cl_Area
## 177 117.2776
## 178 117.6170
## 179 101.9715
## 180 110.1445
## 181 109.4211
## 182 194.9897

# Constants needed for calculations:
clmw <- 35.45      #molecular weight of Chloride: 35.45
smw <- 32.06        #molecular weight of sulfur: 32.06

# Convert ppm to mmol/L
sampledat$S04_mM <- (sampledat$S04_ppm / smw)
sampledat$Cl_mM <- (sampledat$Cl_ppm / clmw)

# Calculate Salinity
# calculated using the Knudsen equation
# Salinity = 0.03 + 1.8050 * Chlorinity
# Ref: A Practical Handbook of Seawater Analysis by Strickland & Parsons (P. 11)

```

```
# =((1.807*Cl_ppm)+0.026)/1000
sampledat$salinity <- ((1.8070 * sampledat$Cl_ppm) + 0.026) / 1000
```

```
head(sampledat)
```

```
##                                     Sample_ID  S04_ppm S04_Area   Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182          TEMPEST_FW_B4_20230606    2.2280  2.0763 173.8469
##   Cl_Area     S04_mM     Cl_mM salinity
## 177 117.2776 17.3441266 117.981489 7.5577019
## 178 117.6170 17.3936619 118.322945 7.5795750
## 179 101.9715 14.9939021 102.583532 6.5713373
## 180 110.1445 16.1900031 110.805585 7.0980268
## 181 109.4211 16.1312757 110.077822 7.0514077
## 182 194.9897  0.0694947   4.904003 0.3141673
```

Pull out dups and check with percent difference

```
#Show me the data that we have from the calculations
head(sampledat)
```

```
##                                     Sample_ID  S04_ppm S04_Area   Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182          TEMPEST_FW_B4_20230606    2.2280  2.0763 173.8469
##   Cl_Area     S04_mM     Cl_mM salinity
## 177 117.2776 17.3441266 117.981489 7.5577019
## 178 117.6170 17.3936619 118.322945 7.5795750
## 179 101.9715 14.9939021 102.583532 6.5713373
## 180 110.1445 16.1900031 110.805585 7.0980268
## 181 109.4211 16.1312757 110.077822 7.0514077
## 182 194.9897  0.0694947   4.904003 0.3141673
```

```
#pull out any rows that have "dup" in the SampleID column
dups <- sampledat %>%
  filter(str_detect(Sample_ID, "_dup"))      #have to change this to match data
head(dups)
```

```
##                                     Sample_ID  S04_ppm S04_Area   Cl_ppm Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.259 109.4211
## 2          TEMPEST_SW_I5_20230606_dup 240.9080  5.6127 2335.150  65.4786
##   S04_mM     Cl_mM salinity
```

```

## 1 16.131276 110.07782 7.051408
## 2 7.514286 65.87165 4.219642

#remove these from sample dataframe in a new dataframe
sampledat2 <- sampledat %>%
  filter(!str_detect(Sample_ID, "_dup")) %>%
  filter(!str_detect(Sample_ID, "_spk"))
head(sampledat2)

##                                     Sample_ID  SO4_ppm SO4_Area   Cl_ppm  Cl_Area
## 1 TEMPEST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438 117.2776
## 2 TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484 117.6170
## 3 TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862 101.9715
## 4 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580 110.1445
## 5 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469 194.9897
## 6 TEMPEST_FW_C3_20230606 12.2001 12.9198 27.3401 34.8467
##           SO4_mM      Cl_mM  salinity
## 1 17.3441266 117.9814894 7.55770195
## 2 17.3936619 118.3229450 7.57957496
## 3 14.9939021 102.5835317 6.57133726
## 4 16.1900031 110.8055853 7.09802681
## 5 0.0694947 4.9040028 0.31416735
## 6 0.3805396 0.7712299 0.04942956

#remove the dup from these IDs so we will have duplicate sample names
dups$Sample_ID<-gsub("_dup","",as.character(dups$Sample_ID))
dups <- dups[, -c(2:5)]
colnames(dups) <- c('Sample_ID', 'SO4_mM_dup', 'Cl_mM_dup', "salinity_dup")
head(dups)

##                                     Sample_ID  SO4_mM_dup Cl_mM_dup salinity_dup
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 16.131276 110.07782 7.051408
## 2 TEMPEST_SW_I5_20230606 7.514286 65.87165 4.219642

#put it back together with the old data set and look for duplicates
QAdups <- merge(sampledat2, dups)
head(QAdups)

##                                     Sample_ID  SO4_ppm SO4_Area   Cl_ppm  Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.058 110.1445
## 2 TEMPEST_SW_I5_20230606 236.2910 5.5051 2285.850 64.0963
##           SO4_mM      Cl_mM  salinity SO4_mM_dup Cl_mM_dup salinity_dup
## 1 16.190003 110.80559 7.098027 16.131276 110.07782 7.051408
## 2 7.370274 64.48097 4.130558 7.514286 65.87165 4.219642

QAdups$SO4_dups_chk <- ((abs(QAdups$SO4_mM-QAdups$SO4_mM_dup))/((QAdups$SO4_mM+QAdups$SO4_mM_dup)/2))*100
QAdups$SO4_dups_flag <- ifelse(QAdups$SO4_dups_chk <10, 'YES', 'NO, rerun')

QAdups$Cl_dups_chk <- ((abs(QAdups$Cl_mM-QAdups$Cl_mM_dup))/((QAdups$Cl_mM+QAdups$Cl_mM_dup)/2))*100
QAdups$Cl_dups_flag <- ifelse(QAdups$Cl_dups_chk <10, 'YES', 'NO, rerun')

head(QAdups)

```

```

##                                     Sample_ID  S04_ppm S04_Area   Cl_ppm  Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.058 110.1445
## 2                 TEMPEST_SW_I5_20230606 236.2910  5.5051 2285.850  64.0963
##   S04_mM      Cl_mM salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk
## 1 16.190003 110.80559 7.098027 16.131276 110.07782    7.051408  0.3633977
## 2 7.370274  64.48097 4.130558  7.514286  65.87165     4.219642  1.9350418
##   S04_dups_flag Cl_dups_chk Cl_dups_flag
## 1          YES  0.6589567      YES
## 2          YES  2.1337241      YES


#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_chk, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call `lifecycle::last_lifecycle_warnings()` to see where this warning was
## generated.

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="Difference Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES        2  YES        2     2    100      100
```

Pull out dups and check with cv

```
#the cv
# calculate the sd bewteen the two duplicates then divide by the mean and multiply by 100
```

```

df2 <- as.data.frame(QAdups$S04_mM)
df2$dups <- QAdups$S04_mM_dup

df2$sds <- apply(df2, 1, sd)

QAdups$S04_dups_cv <- (df2$sds) / ((QAdups$S04_mM + QAdups$S04_mM_dup) / 2) * 100
QAdups$S04_dups_cv_flag <- ifelse(QAdups$S04_dups_cv < 11, 'YES', 'NO, rerun')

head(QAdups)


##                                     Sample_ID  S04_ppm  S04_Area   Cl_ppm  Cl_Area
## 1 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.058 110.1445
## 2                      TEMPEST_SW_I5_20230606 236.2910  5.5051 2285.850  64.0963
##   S04_mM      Cl_mM salinity S04_mM_dup Cl_mM_dup salinity_dup S04_dups_chk
## 1 16.190003 110.80559 7.098027 16.131276 110.07782    7.051408  0.3633977
## 2 7.370274  64.48097 4.130558  7.514286  65.87165    4.219642  1.9350418
##   S04_dups_flag Cl_dups_chk Cl_dups_flag S04_dups_cv S04_dups_cv_flag
## 1          YES     0.6589567        YES     0.256961           YES
## 2          YES     2.1337241        YES     1.368281           YES

#plot dups output as a bar graph to easily check - want any over 10% to be red need to work on this
Sdupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = S04_dups_cv, fill=S04_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between S04 Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

Cldupsbar <- ggplot(data = QAdups, aes(x = Sample_ID, y = Cl_dups_chk, fill=Cl_dups_chk)) +
  geom_bar(stat = 'identity') +
  scale_fill_gradient2(low='darkgreen', mid='darkgreen', high='darkgrey', space='Lab') +
  theme_classic() + labs(x= "Sample ID", y="CV Between Cl Duplicates (%)") +
  theme(legend.position="none") + geom_hline(yintercept=10, linetype="dashed",
                                             color = "black", size=1)

ggarrange(Sdupsbar, Cldupsbar, ncol=2, nrow=1)

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_dups <- QAdups %>%
  group_by(SO4_dups_flag) %>%
  summarise(S_no_rows = length(SO4_dups_flag))
Perc_dups1 <- QAdups %>%
  group_by(Cl_dups_flag) %>%
  summarise(Cl_no_rows = length(Cl_dups_flag))
colnames(Perc_dups) <- c("Flag", "S_no_rows")
colnames(Perc_dups1) <- c("Flag", "Cl_no_rows")
Perc_dups2 <- cbind(Perc_dups, Perc_dups1)

Perc_dups2$Total <- length(QAdups$SO4_dups_flag)
Perc_dups2$S_Percent <- (Perc_dups2$S_no_rows / Perc_dups2$Total)*100
Perc_dups2$Cl_Percent <- (Perc_dups2$Cl_no_rows / Perc_dups2$Total)*100
head(Perc_dups2)
```

```
##   Flag S_no_rows Flag Cl_no_rows Total S_Percent Cl_Percent
## 1  YES        2    YES        2     2      100      100
```

Pull out spikes and check with dionex output conc.

```
#Show me the data that we have from the calculations
head(sampledat)
```

```

##                               Sample_ID  S04_ppm S04_Area   Cl_ppm
## 177      TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438
## 178      TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484
## 179      TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862
## 180      TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182          TEMPEST_FW_B4_20230606    2.2280  2.0763 173.8469
##           Cl_Area     S04_mM     Cl_mM salinity
## 177 117.2776 17.3441266 117.981489 7.5577019
## 178 117.6170 17.3936619 118.322945 7.5795750
## 179 101.9715 14.9939021 102.583532 6.5713373
## 180 110.1445 16.1900031 110.805585 7.0980268
## 181 109.4211 16.1312757 110.077822 7.0514077
## 182 194.9897  0.0694947  4.904003 0.3141673

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk"))      #have to change this to match data
head(spks)

```

```

##                               Sample_ID  S04_ppm S04_Area   Cl_ppm Cl_Area
## 1          TEMPEST_FW_I5_20230606_spk 16.5901 15.4606 31.4901 35.3198
## 2 TEMPEST_FW_SOURCE_20230606_1511_spk 16.3183 17.2810 24.0314 30.6296
##           S04_mM     Cl_mM salinity
## 1 0.5174704 0.8882962 0.05692861
## 2 0.5089925 0.6778956 0.04345074

```

```

#remove the dup from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[,-c(2,3,4,5,7,8)]
colnames(spks) <- c('Sample_ID', 'S04_mM_spk')
head(spks)

```

```

##                               Sample_ID  S04_mM_spk
## 1          TEMPEST_FW_I5_20230606 0.5174704
## 2 TEMPEST_FW_SOURCE_20230606_1511 0.5089925

```

```

#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)

```

```

##                               Sample_ID  S04_ppm S04_Area   Cl_ppm Cl_Area     S04_mM
## 1          TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##           Cl_mM salinity S04_mM_spk
## 1 0.8921834 0.05717762 0.5174704
## 2 0.6922793 0.04437213 0.5089925

```

```

#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250/smw)      # in mM
spkvol <- 10                # in uL
spkvol <- spkvol/1000000

```

```

#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol           # mmoles of S04
head(QAspks)

##                                     Sample_ID S04_ppm S04_Area  Cl_ppm Cl_Area    S04_mM
## 1          TEMPEST_FW_I5_20230606 12.8534  11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474  14.7702 24.5413 31.2795 0.4350405
##      Cl_mM   salinity S04_mM_spk S04_spk_Conc
## 1 0.8921834 0.05717762  0.5174704 7.797879e-05
## 2 0.6922793 0.04437213  0.5089925 7.797879e-05

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
QAspks$Dilution <- 1
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
#QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000000
head(QAspks)

##                                     Sample_ID S04_ppm S04_Area  Cl_ppm Cl_Area    S04_mM
## 1          TEMPEST_FW_I5_20230606 12.8534  11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474  14.7702 24.5413 31.2795 0.4350405
##      Cl_mM   salinity S04_mM_spk S04_spk_Conc Dilution SampleVol
## 1 0.8921834 0.05717762  0.5174704 7.797879e-05           1 1e-06
## 2 0.6922793 0.04437213  0.5089925 7.797879e-05           1 1e-06

```

```

#gives us the total SO4 in the sample in mmoles
QAspks$SO4_Total_unspkd <- (QAspks$SO4_mM/QAspks$Dilution)*(QAspks$SampleVol)

##total SO4 in spiked sample in mmoles
QAspks$SO4_Total_spkd <- (QAspks$SO4_mM_spk/QAspks$Dilution)*(QAspks$SampleVol+spkvol)

QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)
QAspks$spk_recovery <- (QAspks$SO4_Total_spkd/QAspks$SO4_expctd_spkd)*100
QAspks$SO4_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES' , 'NO, rerun'


head(QAspks)

##                                         Sample_ID SO4_ppm SO4_Area   Cl_ppm Cl_Area    SO4_mM
## 1             TEMPEST_FW_I5_20230606 12.8534  11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474  14.7702 24.5413 31.2795 0.4350405
##   Cl_mM   salinity SO4_mM_spk SO4_spk_Conc Dilution SampleVol
## 1 0.8921834 0.05717762  0.5174704 7.797879e-05      1     1e-06
## 2 0.6922793 0.04437213  0.5089925 7.797879e-05      1     1e-06
##   SO4_Total_unspkd SO4_Total_spkd SO4_expctd_spkd spk_recovery SO4_spks_flag
## 1      4.009170e-07   5.692174e-06   7.837971e-05    7.262306    NO, rerun
## 2      4.350405e-07   5.598918e-06   7.841383e-05    7.140217    NO, rerun

#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)

spksbar

```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(SO4_spks_flag) %>%
  summarise(no_rows = length(SO4_spks_flag))
Perc_spks$Total <- length(QAspks$SO4_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   SO4_spks_flag no_rows Total Percent
##   <chr>          <int> <int>   <dbl>
## 1 NO, rerun      2        2     100
```

Pull out spikes and check with area calc

```
#Show me the data that we have from the calculations
head(sampledat)
```

	Sample_ID	SO4_ppm	SO4_Area	Cl_ppm
## 177	TEMPEST_SW_SOURCE_20230606_1301	556.0527	12.9549	4182.4438
## 178	TEMPEST_ESTUARY_SOURCE_20230606_0632	557.6408	12.9919	4194.5484
## 179	TEMPEST_ESTUARY_SOURCE_20230606_1236	480.7045	11.1994	3636.5862
## 180	TEMPEST_ESTUARY_SOURCE_20230606_1524	519.0515	12.0928	3928.0580

```

## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469
## Cl_Area S04_mM Cl_mM salinity
## 177 117.2776 17.3441266 117.981489 7.5577019
## 178 117.6170 17.3936619 118.322945 7.5795750
## 179 101.9715 14.9939021 102.583532 6.5713373
## 180 110.1445 16.1900031 110.805585 7.0980268
## 181 109.4211 16.1312757 110.077822 7.0514077
## 182 194.9897 0.0694947 4.904003 0.3141673

```

```

#now use area of all samples and calculate undilution corrected concentrations in ug/mL
sampledat$S04_ugmL <- ((sampledat$S04_Area-S04_Int)/S04_Slope)
sampledat$Cl_ugmL <- (sampledat$Cl_Area-Cl_Int)/Cl_Slope
head(sampledat)

```

```

## Sample_ID S04_ppm S04_Area Cl_ppm
## 177 TEMEPST_SW_SOURCE_20230606_1301 556.0527 12.9549 4182.4438
## 178 TEMPEST_ESTUARY_SOURCE_20230606_0632 557.6408 12.9919 4194.5484
## 179 TEMPEST_ESTUARY_SOURCE_20230606_1236 480.7045 11.1994 3636.5862
## 180 TEMPEST_ESTUARY_SOURCE_20230606_1524 519.0515 12.0928 3928.0580
## 181 TEMPEST_ESTUARY_SOURCE_20230606_1524_dup 517.1687 12.0490 3902.2588
## 182 TEMPEST_FW_B4_20230606 2.2280 2.0763 173.8469
## Cl_Area S04_mM Cl_mM salinity S04_ugmL Cl_ugmL
## 177 117.2776 17.3441266 117.981489 7.5577019 11.121064 83.64890
## 178 117.6170 17.3936619 118.322945 7.5795750 11.152826 83.89098
## 179 101.9715 14.9939021 102.583532 6.5713373 9.614066 72.73174
## 180 110.1445 16.1900031 110.805585 7.0980268 10.381000 78.56118
## 181 109.4211 16.1312757 110.077822 7.0514077 10.343400 78.04521
## 182 194.9897 0.0694947 4.904003 0.3141673 1.782401 139.07750

```

```

#pull out any rows that have "spk" in the SampleID column
spks <- sampledat %>%
  filter(str_detect(Sample_ID, "spk"))      #have to change this to match data
head(spks)

```

```

## Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area
## 1 TEMPEST_FW_I5_20230606_spk 16.5901 15.4606 31.4901 35.3198
## 2 TEMPEST_FW_SOURCE_20230606_1511_spk 16.3183 17.2810 24.0314 30.6296
## S04_mM Cl_mM salinity S04_ugmL Cl_ugmL
## 1 0.5174704 0.8882962 0.05692861 13.27207 25.19204
## 2 0.5089925 0.6778956 0.04345074 14.83478 21.84672

```

```

## Only keep the columns that we need
spks <- spks[,c(1,9,10)] # dont need this here
head(spks)

```

```

## Sample_ID S04_ugmL Cl_ugmL
## 1 TEMPEST_FW_I5_20230606_spk 13.27207 25.19204
## 2 TEMPEST_FW_SOURCE_20230606_1511_spk 14.83478 21.84672

```

```

#remove the spk from these IDs so we will have duplicate sample names
spks$Sample_ID<-gsub("_spk","",as.character(spks$Sample_ID))
spks <- spks[,c(1,2)]
colnames(spks) <- c('Sample_ID', 'S04_ugmL_spk')
head(spks)

##                                     Sample_ID S04_ugmL_spk
## 1           TEMPEST_FW_I5_20230606    13.27207
## 2 TEMPEST_FW_SOURCE_20230606_1511    14.83478

#put it back together with the old data set and look for duplicates
QAspks <- merge(sampledat, spks)
head(QAspks)

##                                     Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1           TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##          Cl_mM   salinity S04_ugmL Cl_ugmL S04_ugmL_spk
## 1 0.8921834 0.05717762 10.28271 25.30230      13.27207
## 2 0.6922793 0.04437213 12.67940 22.31027      14.83478

#now we need to calculate the spike concentration and calculate the spike recovery
spkconc <- (250)      # in ug
spkvol <- 10            # in uL
spkvol <- spkvol/1000
#spike for these samples was 10uL of the 250mM standard
QAspks$S04_spk_Conc <- (spkconc)*spkvol      # mmoles of S04
head(QAspks)

##                                     Sample_ID S04_ppm S04_Area Cl_ppm Cl_Area   S04_mM
## 1           TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##          Cl_mM   salinity S04_ugmL Cl_ugmL S04_ugmL_spk S04_spk_Conc
## 1 0.8921834 0.05717762 10.28271 25.30230      13.27207      2.5
## 2 0.6922793 0.04437213 12.67940 22.31027      14.83478      2.5

#need to determine dilution factors and initial amount of sample added
#if your samples are all the same dilution just use the first line
#for Steph / COMPASS this depends on the site so...
#QAspks$Dilution <- 1
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 50, QAspks$Dilution)
#$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 50, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 100, QAspks$Dilution)
#$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 50, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 50, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 100, QAspks$Dilution)
#$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 200, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$Dilution)
#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$Dilution)

```

```

#QAspks$Dilution <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$Dilution)

#Set Sample volumes in uL
QAspks$SampleVol <- 1
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "MSM_PPR_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_UP"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_TR"), 1501, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GCrew_WC"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_UP"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_TR"), 1475, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "GWI_WC"), 1462, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_UP"), 100, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_TR"), 100, QAspks$SampleVol)
QAspks$SampleVol <- ifelse(str_detect(QAspks$Sample_ID, "SWH_WC"), 100, QAspks$SampleVol)

#change sample volume to L
QAspks$SampleVol <- QAspks$SampleVol/1000
head(QAspks)

```

```

##                                     Sample_ID SO4_ppm SO4_Area Cl_ppm Cl_Area   SO4_mM
## 1           TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##          Cl_mM  salinity SO4_ugmL Cl_ugmL SO4_ugmL_spk SO4_spk_Conc SampleVol
## 1 0.8921834 0.05717762 10.28271 25.30230      13.27207      2.5    0.001
## 2 0.6922793 0.04437213 12.67940 22.31027      14.83478      2.5    0.001

```

```

#gives us the total SO4 in the sample in mmoles
QAspks$SO4_Total_unspkd <- QAspks$SO4_ugmL*QAspks$SampleVol

##total SO4 in spiked sample in mmoles
QAspks$SO4_Total_spkd <- (QAspks$SO4_ugmL_spk)*(QAspks$SampleVol+spkvol)

QAspks$SO4_expctd_spkd <- (QAspks$SO4_Total_unspkd + QAspks$SO4_spk_Conc)
QAspks$spk_recovery <- (QAspks$SO4_Total_spkd/QAspks$SO4_expctd_spkd)*100
QAspks$SO4_spks_flag <- ifelse(QAspks$spk_recovery <=120 & QAspks$spk_recovery >=80 , 'YES', 'NO, rerun'

head(QAspks)

```


```

##                                     Sample_ID SO4_ppm SO4_Area Cl_ppm Cl_Area   SO4_mM
## 1           TEMPEST_FW_I5_20230606 12.8534 11.9783 31.6279 35.4744 0.4009170
## 2 TEMPEST_FW_SOURCE_20230606_1511 13.9474 14.7702 24.5413 31.2795 0.4350405
##          Cl_mM  salinity SO4_ugmL Cl_ugmL SO4_ugmL_spk SO4_spk_Conc SampleVol
## 1 0.8921834 0.05717762 10.28271 25.30230      13.27207      2.5    0.001
## 2 0.6922793 0.04437213 12.67940 22.31027      14.83478      2.5    0.001
##          SO4_Total_unspkd SO4_Total_spkd SO4_expctd_spkd spk_recovery SO4_spks_flag
## 1          0.01028271      0.1459927      2.510283      5.815788      NO, rerun
## 2          0.01267940      0.1631825      2.512679      6.494364      NO, rerun

```

```
#plot spk recoveries output as a bar graph to easily check - want any over 10% to be red need to work on this
spksbar <- ggplot(data = QAspks, aes(x = Sample_ID, y = spk_recovery)) +
  geom_bar(stat = 'identity') +
  scale_fill_manual(values=c("lightblue4")) +
  theme_classic() + labs(x= "Sample ID", y="Spike Recovery (%)") +
  theme(legend.position="none") +
  geom_hline(yintercept=80, linetype="dashed", color = "black", size=1) +
  geom_hline(yintercept=120, linetype="dashed", color = "black", size=1)

spksbar
```



```
#check for percent of no, reruns to see if it would warrant reruns
Perc_spks <- QAspks %>%
  group_by(S04_spks_flag) %>%
  summarise(no_rows = length(S04_spks_flag))
Perc_spks$Total <- length(QAspks$S04_spks_flag)
Perc_spks$Percent <- (Perc_spks$no_rows / Perc_spks$Total)*100
head(Perc_spks)
```

```
## # A tibble: 1 x 4
##   S04_spks_flag no_rows Total Percent
##   <chr>           <int> <int>   <dbl>
## 1 NO, rerun       2     2     100
```

Make final dataframe with IDs

```
#for steph <- pull out identifiers of the sample names
#pull the sample ID and separate it by the underscores
IDs <- data.frame(do.call('rbind', strsplit(as.character(sampledat2$Sample_ID), '_', fixed=TRUE)))

## Warning in rbind(c("TEMEPST", "SW", "SOURCE", "20230606", "1301"), c("TEMPEST",
## : number of columns of result is not a multiple of vector length (arg 5)

colnames(IDs) <- c("Project", "Plot", "Grid", "Date", "Time")
head(IDs)

##   Project   Plot   Grid   Date   Time
## 1 TEMEPST     SW SOURCE 20230606    1301
## 2 TEMPEST ESTUARY SOURCE 20230606    0632
## 3 TEMPEST ESTUARY SOURCE 20230606    1236
## 4 TEMPEST ESTUARY SOURCE 20230606    1524
## 5 TEMPEST      FW      B4 20230606 TEMPEST
## 6 TEMPEST      FW      C3 20230606 TEMPEST

#rejoin them to the dataframe
alldat <- cbind(IDs, sampledat2)
head(alldat)

##   Project   Plot   Grid   Date   Time           Sample_ID
## 1 TEMEPST     SW SOURCE 20230606    1301 TEMEPST_SW_SOURCE_20230606_1301
## 2 TEMPEST ESTUARY SOURCE 20230606    0632 TEMPEST_ESTUARY_SOURCE_20230606_0632
## 3 TEMPEST ESTUARY SOURCE 20230606    1236 TEMPEST_ESTUARY_SOURCE_20230606_1236
## 4 TEMPEST ESTUARY SOURCE 20230606    1524 TEMPEST_ESTUARY_SOURCE_20230606_1524
## 5 TEMPEST      FW      B4 20230606 TEMPEST          TEMPEST_FW_B4_20230606
## 6 TEMPEST      FW      C3 20230606 TEMPEST          TEMPEST_FW_C3_20230606
##   SO4_ppm SO4_Area Cl_ppm Cl_Area SO4_mM Cl_mM salinity
## 1 556.0527 12.9549 4182.4438 117.2776 17.3441266 117.9814894 7.55770195
## 2 557.6408 12.9919 4194.5484 117.6170 17.3936619 118.3229450 7.57957496
## 3 480.7045 11.1994 3636.5862 101.9715 14.9939021 102.5835317 6.57133726
## 4 519.0515 12.0928 3928.0580 110.1445 16.1900031 110.8055853 7.09802681
## 5  2.2280   2.0763 173.8469 194.9897  0.0694947  4.9040028 0.31416735
## 6 12.2001  12.9198  27.3401  34.8467  0.3805396  0.7712299 0.04942956
```

Make final dataframe with IDs

```
write.csv(alldat, file="Processed Data/COMPASS_TEMPEST_Dionex_ClSo4_Final_20230607.csv") #Char
```

END