COMPASS TEMPEST Discrete DOC Data Workflow: 202505

May 2025

2025-06-05

Run Information

```
#identify which section you are in
cat("Run Information")
```

Run Information

```
#a link to the Gitbook or whatever protocol you are using for this analysis
 #steph will add this soon
#anything that needs to be changed do this in the first chunk
 Date Run = "05/22/25"
 Run_by = "Stephanie J. Wilson"
 Script_run_by = "Stephanie J. Wilson"
 run_notes = "This run, the high check standard was remade halfway
 through run as a 60mg/L C, so checks are okay."
 #file path and name for summary file
   raw_file_name = "tmp_doc_raw_data_2025/TMP_202505.txt"
 #file path and name for the all peaks file
   raw_allpeaks_name = "tmp_doc_raw_data_2025/TMP_202505_allpeaks.txt"
 #file path and name for processed file
   processed_file_name = "tmp_doc_processed_data_2025/TMP_PW_DOC_Processed_202505.csv"
#check standard concentrations - Update if running different checks:
  chk_std_c = 50
  chk std n = 2
#Log path
   Log_path = "tmp_doc_raw_data_2025/COMPASS_TMP_TOCTN_QAQClog_2025.csv"
```

Setup

Pull in active porewater tracking inventory sheet

File already exists. No download needed.

Import Data Functions

Import Sample Data


```
## Import Sample Data
## New names:
## * '' -> '...14'
## # A tibble: 6 x 4
   sample_name
                         npoc_raw tdn_raw run_datetime
##
     <chr>
                          <dbl> <dbl> <chr>
## 1 TMP_CTRL_H3_20250509 25.5 0.977 5/23/2025 2:18:11 AM
## 2 TMP_CTRL_H6_20250509 29.9 0.975 5/23/2025 2:36:46 AM
## 3 TMP_CTRL_I5_20250509 21.5 0.568 5/23/2025 3:06:35 AM
                          10.3 0.365 5/23/2025 3:36:49 AM
18.4 0.657 5/23/2025 3:58:38 AM
## 4 TMP_FW_C6_20250509
## 5 TMP_FW_D5_20250509
                             13.9 0.748 5/23/2025 4:26:17 AM
## 6 TMP_FW_H3_20250509
```

Assessing standard Curves

Assess the Standard Curve


```
## New names:
## 'geom_smooth()' using formula = 'y ~ x'
## * '' -> '...18'
```

NPOC Std Curve by Date




'geom_smooth()' using formula = 'y ~ x'

TN Std Curve by Date



Slope Drift Assessment



- ## [1] "NPOC Curve r2 GOOD"
- ## [1] "TN Curve r2 GOOD"


Assess Check Standards

Assess the Check Standards

```
## New names:
## * '' -> '...14'
```

[1] "Carbon CHECK STANDARD RSD TOO HIGH - REASSESS"

[1] "Nitrogen CHECK STANDARD RSD TOO HIGH - REASSESS"



[1] "<60% of Carbon Check Standards are within range of the expected concentration - REASSESS"

[1] "<60% of Nitrogen Check Standards are within range of the expected concentration - REASSESS"

Assess Blanks


Assess Blanks

New names:

* '' -> '...14'

[1] ">60% of Carbon Blank concentrations are below the lower 25% quartile of samples"

[1] ">60% of Nitrogen Blank concentrations are below the lower 25% quartile of samples"



carbon blanks:

[1] 1.37688


nitrogen blanks:

[1] 0.025593

Assess Duplicates - if there are any

Assess Duplicates


```
## # A tibble: 1 x 3
##
     sample name
                        npoc_raw_dup tdn_raw_dup
     <chr>>
                               <dbl>
                                           <dbl>
## 1 TMP_FW_I5_20250509
                                19.0
                                           0.710
                                                  run_datetime
            sample_name npoc_raw tdn_raw
## 1 TMP_FW_I5_20250509
                            19.3 0.7764 5/23/2025 5:29:16 AM
                                            tdn_flag npoc_raw_dup tdn_raw_dup
                    npoc_flag
## 1 NPOC checks out of range TN checks out of range
                                                             18.96
                                                                        0.7099
            sample_name npoc_raw tdn_raw
                                                 run_datetime
## 1 TMP_FW_I5_20250509
                            19.3 0.7764 5/23/2025 5:29:16 AM
                                            tdn_flag npoc_raw_dup tdn_raw_dup
                    npoc_flag
## 1 NPOC checks out of range TN checks out of range
                                                             18.96
                                                                        0.7099
    npoc_dups_cv npoc_dups_cv_flag tdn_dups_cv tdn_dups_cv_flag
         1.873354
                                YES
                                       9.200138
## 1
## Warning: Using 'size' aesthetic for lines was deprecated in ggplot2 3.4.0.
## i Please use 'linewidth' instead.
## This warning is displayed once every 8 hours.
## Call 'lifecycle::last_lifecycle_warnings()' to see where this warning was
## generated.
```




- ## [1] ">60% of Carbon Duplicates have a CV <10%"
- ## [1] ">60% of Nitrogen Duplicates have a CV <10%"

Sample Flagging

Sample Flagging







Visualize Data by Plot

Visualize Data ## Warning in rbind(c("TMP", "CTRL", "H3", "20250509"), c("TMP", "CTRL", "H6", : ## number of columns of result is not a multiple of vector length (arg 1) Site_Code Plot Grid_Square Date NA TMP CTRL H3 20250509 TMP ## 2 TMP CTRL H6 20250509 TMP TMP CTRL I5 20250509 TMP ## 3 ## 4 TMP FW C6 20250509 TMP ## 5 TMP FW D5 20250509 TMP TMP FW ## 6 H3 20250509 TMP ## Site Code Plot Grid Square Date NA sample name npoc raw tdn raw TMP CTRL H3 20250509 TMP TMP_CTRL_H3_20250509 25.46 0.9769 TMP CTRL H6 20250509 TMP TMP_CTRL_H6_20250509 ## 2 29.91 0.9752 ## 3 TMP CTRL I5 20250509 TMP TMP_CTRL_I5_20250509 21.48 0.5675 ## 4 TMP FW C6 20250509 TMP TMP_FW_C6_20250509 10.30 0.3652 ## 5 TMP FW D5 20250509 TMP TMP_FW_D5_20250509 18.45 0.6574 TMP FW H3 20250509 TMP TMP_FW_H3_20250509 ## 6 13.93 0.7482 run_datetime npoc_flag tdn_flag ## 1 5/23/2025 2:18:11 AM NPOC checks out of range TN checks out of range ## 2 5/23/2025 2:36:46 AM NPOC checks out of range TN checks out of range ## 3 5/23/2025 3:06:35 AM NPOC checks out of range TN checks out of range ## 4 5/23/2025 3:36:49 AM NPOC checks out of range TN checks out of range ## 5 5/23/2025 3:58:38 AM NPOC checks out of range TN checks out of range ## 6 5/23/2025 4:26:17 AM NPOC checks out of range TN checks out of range

Carbon by Plot



Nitrogen by Plot



Convert data from mg/L to uMoles/L

Add in/check metadata

Check Sample IDs with Metadata

```
## # A tibble: 11 x 2
##
      sample_name
                         metadata_recorded
##
      <chr>
                         <1g1>
  1 TMP_C_H3_20250509
##
                        TRUE
## 2 TMP C H6 20250509
                        TRUE
## 3 TMP_C_I5_20250509
                        TRUE
  4 TMP FW C6 20250509 TRUE
## 5 TMP_FW_D5_20250509 TRUE
## 6 TMP_FW_H3_20250509 TRUE
## 7 TMP_FW_H6_20250512 TRUE
## 8 TMP FW I5 20250509 TRUE
## 9 TMP_SW_D5_20250509 TRUE
## 10 TMP_SW_F6_20250509 TRUE
## 11 TMP_SW_H6_20250509 TRUE
```

Export Processed Data

Export Processed Data

```
## # A tibble: 6 x 21
##
                    plot grid Depth_cm sample_type Vial_ID date npoc_mgL npoc_uM
    Project
                    <chr> <chr>
                                   <dbl> <chr>
##
     <chr>>
                                                     <chr>
                                                             <chr>>
                                                                       <dbl>
                                                                               <dbl>
## 1 COMPASS: TEMP~ C
                          НЗ
                                      15 DOC
                                                     C_H3_D~ 2025~
                                                                        25.5
                                                                               2122.
## 2 COMPASS: TEMP~ C
                          Н6
                                      15 DOC
                                                     C H6 D~ 2025~
                                                                        29.9
                                                                               2492.
## 3 COMPASS: TEMP~ C
                          I5
                                      15 DOC
                                                     C I5 D~ 2025~
                                                                        21.5
                                                                               1790
## 4 COMPASS: TEMP~ FW
                          C6
                                      15 DOC
                                                     FW_C6_~ 2025~
                                                                        10.3
                                                                                858.
                                                     FW_D5_~ 2025~
## 5 COMPASS: TEMP~ FW
                          D5
                                      15 DOC
                                                                        18.4
                                                                               1537.
## 6 COMPASS: TEMP~ FW
                          НЗ
                                      15 DOC
                                                     FW_H3_~ 2025~
                                                                        13.9
                                                                               1161.
## # i 12 more variables: npoc_flag <chr>, tdn_mgL <dbl>, tdn_uM <dbl>,
       tdn_flag <chr>, Analysis_runtime <chr>, Run_notes <chr>,
## #
       Evacuation_date_YYYMMDD <dbl>, Collection_Date_YYYYMMDD <dbl>,
       Collection_Start_Time_24hrs <dbl>, Collection_End_Time_24hrs <dbl>,
## #
       EST_EDT <chr>, Volume_mL <dbl>
```

#end