

Tempest: Porewater Nutrients

2025 Samples

2026-05-27

Contents

0.1	Run Information	2
0.2	Assessing standard Curves	3
0.3	Update Standard Log	7
0.4	Dilution Corrections - ensure the latest dilution is kept	9
0.5	Performance Check	9
0.6	Check NOx Reduction Efficiency	10
0.7	Analyze the Check Standards	11
0.8	Analyze Blanks	14
0.9	Analyze Duplicates	18
0.10	Spikes	21
0.11	Matrix Effects	24
0.12	Check to see if samples run match metadata & merge info	28
0.13	Visualize Data by Plot	28

0.1 Run Information

```
cat("Run Information: NAME ") #lets you know what section you're in
```

```
## Run Information: NAME
```

```
#set the run date & user name
sample_year <- "2024|2025|2026"
user <- "Brenden Blakley"

#identify the files you want to read in
#read in as a list to accomodate ultiple runs in a month
NH3_PO4_files<- c("Raw Data/20260504_TMP2025_NP_Run3.csv", "Raw Data/20260508_TMPmix_NP_Run1.csv")
NOx_files      <- c("Raw Data/20260504_TMP2025_NOx_Run3.csv", "Raw Data/20260508_TMPmix_NOx_Run1.csv")

# Define the file path for QAQC log file - NO Need to change just check year
file_path <- "Raw Data/SEAL_TEMPEST_QAQC_Log_2024.csv"
final_path <- "Processed Data/Tempest_Nutrients_2025_Run3_Runmix.csv"

#record any notes about the run or anything other info here:
run_notes <- ""

#Set up file path for metadata
#downloaded metadata csv - downloaded from Google drive as csv for this year
Raw_Metadata = "Raw Data/COMPASS_SynopticCB_PW_SampleLog_2024.csv"

cat(run_notes)
```

```
##Setup
```


```
##Pull in active porewater tracking inventory sheet from Google Drive:
```

```
##Create similar sample IDs to match with run samples
```

```
##Import Data & Clean
```


0.2 Assessing standard Curves

```
## `geom_smooth()` using formula = 'y ~ x'
```




```
## `geom_smooth()` using formula = 'y ~ x'
```

PO4 Standard Curve






NOx Standard Curve



Test	Run_Date	R2_Flag	R2	Slope
Ammonia 2	05/04/2026	R2_Pass	0.9995982	1.818107
Ammonia 2	05/08/2026	R2_Pass	0.9999344	1.758303
o-PHOS 0.3	05/04/2026	R2_Pass	0.9993257	2.405344
o-PHOS 0.3	05/08/2026	R2_Pass	0.9993370	2.335316
Vanadium NOx	05/04/2026	R2_Pass	0.9997660	1.159949
Vanadium NOx	05/08/2026	R2_Pass	0.9989603	1.272420

0.3 Update Standard Log





0.4 Dilution Corrections - ensure the latest dilution is kept

```
## [1] "Reran Samples: TMP_SW_D5_20251101"
```


```
## [1] "Diluted Samples: TMP_SW_D5_20251101"
```

```
## [1] "No Naming Issues Detected"
```

0.5 Performance Check

Test	Run_Date	PE_Flag	PE_Conc	PE_Target_Conc
Ammonia 2	05/04/2026	Performance Check Within 25% - PROCEED	1.2530565	1.034
Ammonia 2	05/08/2026	Performance Check Within 25% - PROCEED	1.2651620	1.034
o-PHOS 0.3	05/04/2026	Performance Check Within 25% - PROCEED	0.9222755	0.824
o-PHOS 0.3	05/08/2026	Performance Check Within 25% - PROCEED	0.9105505	0.824
Vanadium NOx	05/04/2026	Performance Check Within 25% - PROCEED	1.5966770	1.510
Vanadium NOx	05/08/2026	Performance Check Within 25% - PROCEED	1.6844615	1.510


0.6 Check NOx Reduction Efficiency




Test	Run_Date	Red_Eff_Flag
	05/04/2026	Mean NOx Reduction Efficiency >85% - PROCEED
	05/08/2026	Mean NOx Reduction Efficiency >85% - PROCEED

0.7 Analyze the Check Standards


Test	Run_Date	RSV_Flag	RSV	RSV_Cutoff
Ammonia 2	05/04/2026	RSV WITHIN RANGE - PROCEED	0.0207134	0.25
Ammonia 2	05/08/2026	RSV WITHIN RANGE - PROCEED	0.0207134	0.25
o-PHOS 0.3	05/04/2026	RSV WITHIN RANGE - PROCEED	0.0150140	0.25
o-PHOS 0.3	05/08/2026	RSV WITHIN RANGE - PROCEED	0.0150140	0.25
Vanadium NOx	05/04/2026	RSV WITHIN RANGE - PROCEED	0.0269662	0.25
Vanadium NOx	05/08/2026	RSV WITHIN RANGE - PROCEED	0.0269662	0.25



PO4 Check Standards




NOx Check Standards




Test	Run_Date	CHK_Flag
Ammonia 2	05/04/2026	>60% of Checks Pass - PROCEED
Ammonia 2	05/08/2026	>60% of Checks Pass - PROCEED
o-PHOS 0.3	05/04/2026	>60% of Checks Pass - PROCEED
o-PHOS 0.3	05/08/2026	<60% of Checks Pass - REASSESs
Vanadium NOx	05/04/2026	>60% of Checks Pass - PROCEED
Vanadium NOx	05/08/2026	>60% of Checks Pass - PROCEED

0.8 Analyze Blanks


Test	Run_Date	BLK_Pct_Flag	Mean_Blkc_Conc	Quantile_25
Ammonia 2	05/04/2026	>60% of Blanks Pass - PROCEED	-0.0249550	0.0246705
Ammonia 2	05/08/2026	>60% of Blanks Pass - PROCEED	-0.0027873	0.0175657
o-PHOS 0.3	05/04/2026	>60% of Blanks Pass - PROCEED	0.0008352	0.0046075
o-PHOS 0.3	05/08/2026	>60% of Blanks Pass - PROCEED	0.0014364	0.0105687
Vanadium NOx	05/04/2026	>60% of Blanks Pass - PROCEED	-0.0004657	0.0125000
Vanadium NOx	05/08/2026	>60% of Blanks Pass - PROCEED	-0.0022839	0.0125000





PO4 Blanks




NOx Blanks



0.9 Analyze Duplicates






NOx Duplicate Check



Test	Run_Date	Dup_Flags
Ammonia 2	05/04/2026	<60% of Dups Pass - REASSESS
Ammonia 2	05/08/2026	<60% of Dups Pass - REASSESS
o-PHOS 0.3	05/04/2026	>60% of Dups Pass - PROCEED
o-PHOS 0.3	05/08/2026	>60% of Dups Pass - PROCEED
Vanadium NOx	05/04/2026	<60% of Dups Pass - REASSESS
Vanadium NOx	05/08/2026	<60% of Dups Pass - REASSESS



0.10 Spikes






Test	Run_Date	Spike_Flags
Ammonia 2	05/04/2026	>60% of Spikes have a Diff <50% - PROCEED
Ammonia 2	05/08/2026	>60% of Spikes have a Diff <50% - PROCEED
Vanadium NOx	05/04/2026	>60% of Spikes have a Diff <50% - PROCEED
Vanadium NOx	05/08/2026	>60% of Spikes have a Diff <50% - PROCEED
o-PHOS 0.3	05/04/2026	>60% of Spikes have a Diff <50% - PROCEED
o-PHOS 0.3	05/08/2026	>60% of Spikes have a Diff <50% - PROCEED


0.11 Matrix Effects







Test	Run_Date	Matrix_Flags
Ammonia 2	05/04/2026	Matrix Has Diff <50% - PROCEED
Ammonia 2	05/08/2026	Matrix Has Diff <50% - PROCEED
Vanadium NOx	05/04/2026	Matrix Has Diff <50% - PROCEED
Vanadium NOx	05/08/2026	Matrix Has Diff <50% - PROCEED
o-PHOS 0.3	05/04/2026	Matrix Has Diff <50% - PROCEED
o-PHOS 0.3	05/08/2026	Matrix Has Diff <50% - PROCEED

##Sample Flagging - Within range of standard curve







```
##Add QAQC Flags
```

0.12 Check to see if samples run match metadata & merge info

```
## Some sample IDs are missing from metadata.
```

```
## [1] "TMP_C_H6_20260602"
```


```
## [1] "Sample_Name" "Run_Number" "Conc" "Absorbance" "Dilution"
## [6] "Unit" "Test" "Run_Time" "Run_Date" "Keep"
## [11] "Pair_ID" "Conc_uM" "Conc_flag"
```


```
##Format Data
```


0.13 Visualize Data by Plot

```
## Site_Code Plot Grid_Square Date
## 1 TMP C B4 20260309
## 2 TMP C C3 20260309
## 3 TMP C D5 20260309
## 4 TMP C F6 20250711
## 5 TMP C H3 20250711
## 6 TMP C H3 20260309
```

##	Site_Code	Plot	Grid_Square	Date	Sample_ID	Site	Zone	Grid
## 1	TMP	C	B4	20260309	TMP_C_B4_20260309	Tempest	C	B4
## 2	TMP	C	C3	20260309	TMP_C_C3_20260309	Tempest	C	C3
## 3	TMP	C	D5	20260309	TMP_C_D5_20260309	Tempest	C	D5
## 4	TMP	C	F6	20250711	TMP_C_F6_20250711	Tempest	C	F6
## 5	TMP	C	H3	20250711	TMP_C_H3_20250711	Tempest	C	H3
## 6	TMP	C	H3	20260309	TMP_C_H3_20260309	Tempest	C	H3
##	Depth_cm	Event_Time	Time_of_day	Source	Volume_mL	Collection_Date		
## 1	15	<NA>	<NA>	PW	9.0	20260309		
## 2	15	<NA>	<NA>	PW	4.0	20260309		
## 3	15	<NA>	<NA>	PW	5.0	20260309		
## 4	15	<NA>	<NA>	PW	7.5	20250711		
## 5	15	<NA>	<NA>	PW	14.0	20250711		
## 6	15	<NA>	<NA>	PW	8.0	20260309		
##	Evacuation_Date	NOx_Conc_mgL	NOx_Conc_uM	NOx_Conc_Flag	NOx_QAQC_Flag			
## 1	20260306	0.002901	0.2071152		bd1			
## 2	20260306	0.001754	0.1252258		bd1			
## 3	20260306	0.012405	0.8856476		bd1			
## 4	20250708	0.007791	0.5562338		bd1			
## 5	20250708	0.002212	0.1579244		bd1			
## 6	20260306	0.024258	1.7318855		bd1			
##	NH3_Conc_mgL	NH3_Conc_uM	NH3_Conc_Flag	NH3_QAQC_Flag	P04_Conc_mgL	P04_Conc_uM		
## 1	0.028899	2.063227	Within_Range		0.004516	0.1458008		
## 2	0.051426	3.671529	Within_Range		0.018202	0.5876587		
## 3	0.037136	2.651303	Within_Range		0.007024	0.2267726		
## 4	0.065259	4.659127	Within_Range		0.029424	0.9499653		
## 5	0.077177	5.510006	Within_Range		0.013632	0.4401144		
## 6	0.255198	18.219709	Within_Range		0.019172	0.6189755		
##	P04_Conc_Flag	P04_QAQC_Flag	Analysis_rundate	Run_notes	Field_notes			
## 1	Within_Range		05/08/2026		<NA>			
## 2	Within_Range		05/08/2026		<NA>			
## 3	Within_Range		05/08/2026		<NA>			
## 4	Within_Range		05/04/2026		<NA>			
## 5	Within_Range		05/04/2026		<NA>			
## 6	Within_Range		05/08/2026		<NA>			







##Write Out Data

#end